

MedGyn Aspiration Kit Double Valve

The MedGyn Aspiration Kit Double Valve is used for uterine aspiration/uterine evacuation in obstetrics and gynecologic procedures.

Features & Benefits

- For uterine aspiration and endometrial biopsy
- Volume: 60 mL
- Vacuum: 610-660mm (24-26in) Hg
- Made of latex-free plastic
- Individually-packaged and sterilized
- Fits MedGyn 4, 5, 6, 7, 8, 9, 10, 11, 12mm flexible Karman cannula
- Fits a wide range of MedGyn rigid straight and rigid curved curettes



022526 MedGyn Aspiration Kit Double Valve



Flexible Cannula/Karman

- 022004K Flexible Cannula/Karman 4mm
- 022005K Flexible Cannula/Karman 5mm
- 022006K Flexible Cannula/Karman 6mm
- 022007K Flexible Cannula/Karman 7mm

- 022008K Flexible Cannula/Karman 8mm
- 022009K Flexible Cannula/Karman 9mm
- 022010K Flexible Cannula/Karman 10mm
- 022011K Flexible Cannula/Karman 11mm
- 022012 Flexible Cannula/Karman 12mm



Rigid Curved Curette

- 022106 Curette Rigid Curved 6mm
- 022107 Curette Rigid Curved 7mm
- 022108 Curette Rigid Curved 8mm
- 022109 Curette Rigid Curved 9mm

- 022110 Curette Rigid Curved 10mm
- 022111 Curette Rigid Curved 11mm
- 022112 Curette Rigid Curved 12mm



Rigid Straight Curette

- 022206 Curette Rigid Straight 6mm
- 022207 Curette Rigid Straight 7mm
- 022208 Curette Rigid Straight 8mm
- 022209 Curette Rigid Straight 9mm

- 022210 Curette Rigid Straight 10mm
- 022211 Curette Rigid Straight 11mm
- 022212 Curette Rigid Straight 12mm

TRUSTED PRODUCTS FOR PREMIER WOMEN'S HEALTHCARE

MEDGYN PRODUCTS, INC.

100 W. Industrial Rd. Addison, IL 60101 USA
 t: +1 630.627.4105 toll-free: 800.451.9667
 f: +1 630.627.0127
 e: info@medgyn.com
 w: medgyn.com
 ASPIRDOB-PF2

