


Textured Silicone Implants

ROUND BASE

2) Round Profile / Medial Pole / High Projection



Ref	Vol	K	
	cm ³	cm	
		A	B
20621-115	115	8.1	3.2
20621-135	135	8.7	3.3
20621-155	155	8.8	3.6
20621-175	175	9.3	3.7
20621-195	195	9.7	3.9
20621-215	215	10.0	3.9
20621-235	235	10.5	3.9
20621-255	255	11.0	3.9
20621-285	285	11.4	4.0
20621-305	305	11.7	4.0
20621-325	325	11.8	4.4
20621-355	355	12.1	4.4
20621-385	385	12.5	4.5
20621-435	435	13.2	4.6
20621-485	485	13.6	4.8
20621-525	525	13.8	5.0
20621-575	575	14.1	5.2
20621-625	625	14.6	5.3
20621-695	695	15.1	5.5
20621-755	755	15.6	5.9
20621-815	815	16.0	6.0
20621-880	880	16.4	6.1
20621-935	935	16.8	6.3
20621-995	995	17.1	6.4

Autoclavable breast sizers available in all sizes