


Textured Silicone Implants

OVAL BASE

7) Anatomical Profile / Moderate Height / Moderate Projection



Ref	Vol	K		
	cm ³	cm		
		A	B	C
20645-120L	120	8.9	8.0	3.4
20645-150L	150	9.4	8.5	3.6
20645-180L	180	10.0	8.9	3.9
20645-210L	210	11.0	9.6	4.0
20645-250L	250	12.0	10.3	4.2
20645-290L	290	12.3	10.7	4.5
20645-330L	330	13.0	11.2	4.5
20645-370L	370	13.4	11.4	4.7
20645-400L	400	13.7	11.8	4.9
20645-450L	450	14.1	12.2	5.1
20645-500L	500	14.7	12.8	5.3
20645-550L	550	15.2	13.1	5.4

Autoclavable breast sizes available in all sizes