


Gluteal Implants


OVAL BASE

1) Low Projection Opaque Surface



Ref	Vol	K		
	cm ³	cm		
		A	B	C
13960-200	200	14.0	9.5	2.7
13960-250	250	14.5	9.8	3.1
13960-300	300	14.9	10.4	3.5
13960-350	350	15.2	10.9	3.8
13960-400	400	15.4	11.2	4.3
13960-450	450	15.8	11.7	4.5
13960-500	500	16.2	12.1	4.6

2) High Projection Opaque Surface



Ref	Vol	K		
	cm ³	cm		
		A	B	C
13961-200	200	13.7	9.2	2.8
13961-250	250	14.0	9.5	3.5
13961-300	300	14.3	9.7	3.8
13961-350	350	14.9	10.4	4.2
13961-400	400	15.2	10.9	4.6
13961-450	450	15.4	11.2	4.7
13961-500	500	15.8	11.7	4.8

Gluteal Implants

ROUND BASE

1) Single Projection Smooth Surface



—A—|B|

Ref	Vol	K	
	cm ³	cm	
		A	B
10910-200	200	10.6	3.7
10910-220	220	10.7	3.8
10910-240	240	10.9	3.9
10910-270	270	11.4	3.9
10910-300	300	11.7	4.0
10910-330	330	12.1	4.1
10910-360	360	12.3	2.8