



## Rectangle Gewebe-Expander

Order No.	Item	Before Swelling			After Swelling in 0,9% NaCl-Sol.		
		Volume in ml	Projection in mm	Length x Width mm	Volume ml	in Projection in mm	Length x Width in mm
352-2060	Rectangle 60 ml	5	13	25 x 13	60	36	60 x 36
352-2075	Rectangle 75 ml	6	12	32 x 16	75	35	74 x 41
352-2130	Rectangle 130 ml	13	15	40 x 20	130	45	85 x 50
352-2200	Rectangle 200 ml	20	18	45 x 24	200	52	95 x 60
352-2300	Rectangle 300 ml	30	21	54 x 28	300	58	115 x 65



## Round Breast Tissue Expander

Order No.	Item	Before swelling			After swelling in 0,9% NaCl-Sol.		
		Volume ml	in Projection in mm	Diameter in mm	Volume ml	in Projection in mm	Diameter in mm
352-1200	Round 200 ml	20	18	37	200	52	80
352-1330	Round 330 ml	30	21	42	330	49	101
352-1450	Round 450 ml	43	24	51	450	60	110
352-1550	Round 550 ml	60	26	51	550	71	115
352-1650	Round 650 ml	70	28	55	650	75	120



CE 0297

osmed gmbh  
 Ehrenbergstraße 11  
 98693 Ilmenau  
 Fon: 03677/668 631  
[info@osmed.biz](mailto:info@osmed.biz)  
[www.osmed.biz](http://www.osmed.biz)

Durchmesser  
 PA25.08.06