



Blade for Humeca D80 dermatome

Humeca supplies a range of high quality blades for different types of dermatomes.


For the Humeca® D42 and D80 dermatomes and for the SOBER dermatome Humeca supplies symmetrical and double facet grinded blades.

Humeca also supplies blades that fit in Aesculap® / B.Braun® cordless dermatomes, Integra® / Padgett® dermatomes type B, C and S and Zimmer® dermatomes type 8801 (air) and 8821 (electric).

Blades: features

- Precision facet grinded without burrs.
- Material: stainless steel 1.4034.
- Separately packed in pouch with safety sleeve.
- Radiation sterilized.
- Supplied in boxes of 10 pcs.

Please note that these blades are manufactured by Humeca and suppliers of the original brand blades are not responsible for their performance



Humeca blade that fits Aesculap® B. Braun® dermatome



Humeca blade that fits Integra® Padgett® dermatome



Humeca blade that fits Zimmer® dermatome

Ordering

5.D42BL10	D42 dermatome blades, box 10 pcs.
5.D80BL10	D80 dermatome blades, box 10 pcs.
5.BLZM10	Blades for Zimmer Electric and air dermatomes. Box of 10
5.BL4245	Blades for Aesculap® cordless dermatomes, box 10 pcs.
5.BLPG10	Blades for Padgett® dermatomes, model B, C and S, box 10 pcs.
5.BLSB10	Blades for Sober hand dermatome, box 10 pcs.



Humeca introduces a new range of grooved skin graft carriers for expansion and perforation, called ‘V-carriers’.



V - C A R R I E R S F E A T U R E S

- Compatible with existing Zimmer® and Aesculap® / B.Braun® meshers
- Symmetric V-pattern of grooves prevents sideward movement
- Standard length 280 mm (11.0"), width 78.8 mm (3.1")
- Expansion ratio 1:1.5, 1:2 (type V10 only) and 1:3
- Special carrier for 1:1 perforation (drainage) of sheet grafts
- Highly flexible medical grade polypropylene material
- Individually sterile packed in peel pouch; standard boxes 10 pcs.
- Attractively priced

H U M E C A[®] V - C A R R I E R S

6.V10-1.0	V-carrier, perforation 1:1, for Humeca and Zimmer® meshers, box 10 pcs.
6.V10-1.5	V-carrier, expansion ratio 1:1.5, for Humeca and Zimmer® meshers, box 10 pcs.
6.V10-2.0	V-carrier, expansion ratio 1:2, for Humeca and Zimmer® meshers, box 10 pcs.
6.V10-3.0	V-carrier, expansion ratio 1:3, for Humeca and Zimmer® meshers, box 10 pcs.
6.V15-1.0	V-carrier, perforation 1:1, for Humeca and Aesculap® / B.Braun® meshers, box 10 pcs.
6.V15-1.5	V-carrier, expansion ratio 1:1.5, for Humeca and Aesculap® / B.Braun® meshers, box 10 pcs.
6.V15-3.0	V-carrier, expansion ratio 1:3, for Humeca and Aesculap® / B.Braun® meshers, box 10 pcs.