



# Body Contouring & Lipofilling Catalogue

[www.eurosurge.co.uk](http://www.eurosurge.co.uk)



## Body Contouring & Lipofilling Catalogue

Vibrasat®, Liposat® & Vacusat®	3
Consumables	4
Cannula Portfolio for Vibrasat® Power	5
Cannula Portfolio for Vibrasat® Pro	6
Liposuction Equipment - Aspiration	7
Liposuction Cannulae - Single Use	8
Liposuction Cannulae - Re-useable	9
Liposuction Accessories & Additional Items	10 - 11
Infiltration Pumps & Accessories	12 - 13
Macro & Micro Fat Grafting Kits - Single Use	14
AquaVage - High Volume Fat Grafting - Single Use	15
Cannulae For Fat Harvesting - Single Use	16
Cannulae For Injection/Filling - Single Use	17 - 18
Fat Injection Cannulae - Re-usable	19
Female Garments	20
Male Garments	21
Garment Accessories	21
Our Garment Range	22



## Vibrasat® power The Original Vibrasat®

The Vibrasat® power device moves the cannula in very fast and short longitudinal movements. Combined with the suction vacuum, these movements lead to a gentle removal of adipose tissue. The cannula glides through the body with low effort, which reduces strain on the surgeon and provides flexibility for surgeons to easily maneuver the cannula.

Order Ref; 00002246 - Console plus 1 hand piece

Order Ref; 00002327 - Hand piece for Vibrasat power



## Vibrasat® Pro New Vibrating Handle and Console for Liposuction

Vibrasat® Pro console is extremely quiet and reciprocates up to 5000 strokes per minute (SPM) with a 3mm stroke with a boost function of 6000 SPM. New in the Vibrasat® Pro is the QuickLock® connection which makes changing cannulae easy and fast. Also the new button on the handle provides quick access to the boost functions.

Order Ref; 00003920 - Vibrasat® Pro Console plus 1 hand piece

Order Ref; 00003922 - Hand piece for Vibrasat® Pro



## Liposat® power Infiltration Device

Liposuction using the tumescent local anesthesia (TLA) technique is the most frequently performed procedure around the world. The infiltration pump Liposat® power has been specially developed for this application. Its features and functions have been designed and adjusted for surgical and aesthetic medicine. This method has principal advantages for any treatment that requires local anaesthesia, such as liposuction or vascular surgery. Liposat Infiltration pump can be used in conjunction with Aquavage, Soft Touch irrigation tubing and Vibrasat for a complete closed and Power Assisted Fat grafting / body contouring procedure (as described by Del Vecchio & Wall).

Order Ref; 00002274



## Vacusat® power Surgical Aspirator

The Vacusat power is a high-performance surgical aspirator that is precisely designed to meet the needs of modern liposuction. The silent and powerful Vacusat® power fulfills all the specific requirements in the operating room. Thanks to its innovative sound insulation technology, it is one of the quietest suction devices in the market that is easy to operate.

Order Ref; 00002252

## Consumables



**Tumescence Tube  
for Liposat® power**  
Order Ref; 00002251



**Suction Tube (red)**  
Order Ref; 00002264



**Bacteria Filter**  
Order Ref; 00002251  
077.0571



**Hydrophobic Filter**  
Order Ref; 00002297



**Overflow Protection**  
Order Ref; 00002299



**Medela Disposable  
Suction Liners**

Order Ref;

077.0083	1,5 litre suction liner	(box of 20)
077.0084	1,5 litre suction liner with solidifier	(box of 15)
077.0086	2,5 litre suction liner	(box of 20)
077.0087	2,5 litre suction liner with solidifier	(box of 15)
077.0571	Bacteria filter	
077.0563	Medela changeover valve with vacuum tubing	

## Accessories



**Heating Dish**  
Order Ref; 00002542



**Holding Clamp**  
Order Ref; 00002258















**Attachment Kit Vibrasat®**  
Order Ref;  
00002283 (Power) 00003973 (Pro)



**3-Pedal Footswitch**  
Order Ref;  
93003517 (2m) 00003981 (5m)













## Cannulae For Vibrasat Power

Single use cannulae come in boxes of 10. Reusable cannulae are sold individually.

Description	Single Use Order Ref;	Reuseable Order Ref;	Hole	Ø (mm)	Length (mm)
 Single-Port	ASC-S210-30-0350	ASC-R210-30-0350	1	3	350
 Tri-Port	ASC-S230-40-0350	ASC-R230-40-0350	3	4	350
 Multi-Hole	-	ASC-R310-30-0250	24	3	250
	ASC-S310-30-0350	ASC-R310-30-0350	24	3	350
 Mercedes	-	ASC-R010-25-0150	3	2.5	150
	00002859	-	3	2	250
	00002857	ASC-R010-30-0350	3	3	350
	00002858	ASC-R010-30-0250	3	3	250
	00002855	ASC-R010-40-0250	3	4	250
	00002856	ASC-R010-40-0350	3	4	350
 Double Mercedes	ASC-S020-40-0350	ASC-R020-40-0350	6	4	350
 Triple Mercedes	ASC-S030-40-0350	ASC-R030-40-0350	9	4	350
 Rapid Extraction	ASC-S100-25-0150	ASC-R100-25-0150	9	2	150
 9-Hole-Rapid Extraction	00002891	ASC-R100-30-0250	9	3	250
	00002893	ASC-R100-40-0250	9	4	250
	00002892	ASC-R100-40-0350	9	4	350
	-	ASC-R100-30-0350	9	3	350
 Multi Rapid Extraction	ASC-S110-40-0350	ASC-R110-40-0350	18	4	350
 Sweat gland suction	00002889	-	9	3	-
 Basket Cannula	ASC-S510-40-0350	ASC-R510-40-0350	-	4	350
	ASC-S510-50-0350	ASC-R510-50-0350	-	5	350
 Double Basket	ASC-S530-40-0350	ASC-R530-40-0350	-	4	350
	ASC-S530-50-0350	ASC-R530-50-0350	-	5	350

## Cannulae For Vibrasat Pro

Single use cannulae come in boxes of 10. Reusable cannulae are sold individually.

Description	Single Use	Reuseable	Hole	Ø (mm)	Length (mm)
 Basket Cannula	ASC-S511-40-0350	ASC-R511-40-0350		4	350
	ASC-S511-50-0350	ASC-R511-50-0350		5	350
 Candy Cane	ASC-S411-40-0250	ASC-R411-40-0250		4	250
	ASC-S411-40-0350	ASC-R411-40-0350		4	350
 Double Basket	ASC-S531-40-0350	ASC-R531-40-0350		4	350
	ASC-S531-50-0350	ASC-R531-50-0350		5	350
 Double Mercedes	-	ASC-R021-40-0250	6	4	250
	ASC-S021-40-0350	ASC-R021-40-0350	6	4	350
 Mercedes	ASC-S011-25-0150	ASC-R011-25-0150	3	2.5	150
	ASC-S011-30-0250	ASC-R011-30-0250	3	3	250
	ASC-S011-30-0350	ASC-R011-30-0350	3	3	350
	ASC-S011-40-0250	ASC-R011-40-0250	3	4	250
	ASC-S011-40-0350	ASC-R011-40-0350	3	4	350
 Multi Rapid Extraction	ASC-S111-40-0350	ASC-R111-40-0350	18	4	350
	ASC-S311-30-0250	ASC-R311-30-0250	24	3	250
 Multi-Hole	ASC-S311-30-0350	ASC-R311-30-0350	24	3	350
	ASC-S311-40-0350	ASC-R311-40-0350	24	4	350
 Rapid Extraction	ASC-S101-25-0150	ASC-R101-25-0150	9	2.5	150
	ASC-S101-30-0250	ASC-R101-30-0250	9	3	250
	ASC-S101-40-0250	ASC-R101-40-0250	9	4	250
	ASC-S101-40-0350	ASC-R101-40-0350	9	4	350
 Single-Port	ASC-S211-30-0250	ASC-R211-30-0250	1	3	250
	ASC-S211-30-0350	ASC-R211-30-0350	1	3	350
 Tri-Port	ASC-S231-30-0250	ASC-R231-30-0250	3	3	250
	-	ASC-R231-30-0350	3	3	350
	ASC-S231-40-0250	ASC-R231-40-0250	3	4	250
	ASC-S231-40-0350	ASC-R231-40-0350	3	4	350
 Triple Mercedes	ASC-S031-40-0350	ASC-R031-40-0350	9	4	350
 9 Hole Rapid Extraction	ASC-S101-30-0350	ASC-R101-30-0350	9	3	350

## Aspiration Pump

### CA.MI Hospivac 400 Economical Liposuction System

Electrically driven with a high pressure motor, this machine delivers fast and quiet adjustable pressure up to 760mmHg.

Two 3 litre suction containers with change-over valve and foot pedal.

Order Ref; RE410350-25 Hospivac 400



## Accessories

### FloVac Suction Jars & Liners – CA.MI Hospivac

FloVac liners are disposable suction liner bags with rigid lids that fit directly onto the FloVac suction jars. 2 and 3 litre liners are available and are fitted with overflow filters to prevent damage to the suction systems.

Order Ref; 31846 2 litre liner

Order Ref; 31847 3 litre liner

Order Ref; 31841 2 litre jar

Order Ref; 31842 3 litre jar

Order Ref; SP0047 bacterial filter



### Medela - Dominant Flex Suction Pump

One-touch adaptable flow rates allow for a wider variety of procedures to be carried out with just one pump. Designed for use in hospitals, clinics and medical practices, the Dominant Flex is the optimal medical suction pump for all suctioning needs. Fitted with 2 x 2,5 litre bottles and a changeover valve. Optional on/off foot switch is available.

Order Ref; 071.003T

077.0571 Bacteria filter

077.0563 The Medela changeover valve with vacuum tubing

### Medela Disposable Suction Liners

For use with ALL versions of Medela Aspiration pumps, such as Vario, Basic 30 and Dominant 50.

Order Ref;

077.0083 1,5 litre suction liner (box of 20)

077.0084 1,5 litre suction liner with solidifier (box of 15)

077.0086 2,5 litre suction liner (box of 20)

077.0087 2,5 litre suction liner with solidifier (box of 15)



### Suction Tubing

Easy to fit with two flexible connectors that will accommodate all types of push-fit cannulae and connections. Available in clear high strength plastic that prevents collapsing even at high aspiration pressure. Packaged 10 per box.

Order Ref;

ST30200 – 10 200cm tubing ø10mm (external) 6.2mm (internal) with 11.2mm connector (internal)







ST 30300 -10 300cm tubing ø10mm (external) 6.2mm (internal) with 11.2mm connector (internal)

HPF 1780 400cm reinforced corrugated tubing ø 9mm (external) 5.7mm (internal) with 10mm connector (internal)



## Single Use

Single-use, sterile range of Liposuction cannulae with easy to hold pre-fixed handle.  
Available in a wide range of tip styles, length and width. 6 cannulae per box.

		REF.	Hole	Ø (mm)	Length (mm)
MERCEDES		LC31302	2	3,00	130
		LC32302		3,00	230
		LC42302		4,00	230
		LC42802		4,00	280
		LC52802		5,00	280
		LC62802		6,00	280
		LC31303	3	3,00	130
		LC32303		3,00	230
		LC42303		4,00	230
		LC42803		4,00	280
		LC52803		5,00	280
		LC62803		6,00	280
SUB MENTAL		LC31503HVS	3	3,00	150
		LC32303HVS		3,00	230
		LC31304H	4	3,00	130
		LC32804H		3,00	280
FAT HARVEST		LC313024H	24	3,00	130
		LC328024H		3,00	280
		LC428024H		4,00	280
		DOLF3150	3	3,00	150
		DOLF3250		3,00	250
		DOLF4150		4,00	150
		DOLF4250		4,00	250



B 89ASH-LL

Luer Lock Hub

B 89TAH-60CC

60cc Toomey Syringe Hub



Luer Lock Aspiration Handle

Toomey Aspiration "Delrin" Handle

## Tip Styles



Mercedes Pattern

Accelerator Pattern

Vegas Pattern

Standard Two Pattern

Standard Pattern

Candy Cane Pattern

Man Pattern

Lilis Keel Pattern



Spatula Pattern

Fournier Pattern

Basket Pattern

Keel Cobra Pattern

Sattler Tip

Fat Disruptor Tip

Cobra Style Tip

## Vitruvian Liposuction Cannula Popular Patterns

### Mercedes

B 8MER-215L	2 x 150mm	Vitruvian Handle
B 8MER-226L	2 x 260mm	Vitruvian Handle
B 8MER-215-SYR	2 x 150mm	Toomey Hub
B 8MER-326L	3 x 260mm	Vitruvian Handle
B 8MER-326-SYR	3 x 260mm	Toomey Hub
B 8MER-332L	3 x 320mm	Vitruvian Handle
B 8MER-415L	4 x 150mm	Vitruvian Handle
B 8MER-426-SYR	4 x 260mm	Toomey Hub
B 8MER-432L	4 x 320mm	Vitruvian Handle
B 8MER-526L	5 x 260mm	Vitruvian Handle

### Accelerator

B 8III-226L	2 x 260mm	Vitruvian Handle
B 8III-215-SYR	2 x 150mm	Toomey Hub
B 8III-315L	3 x 150mm	Vitruvian Handle
B 8III-326L	3 x 260mm	Vitruvian Handle
B 8III-332L	3 x 320mm	Vitruvian Handle
B 8III-426-SYR	4 x 260mm	Toomey Hub
B 8III-432L	4 x 320mm	Vitruvian Handle

### Fat Disruptor

B 8MBAT-315L	3 x 150mm	Vitruvian Handle
B 8MBAT-326-SYR	3 x 260mm	Toomey Hub

### Cobra

B 8COB-326L	3 x 260mm	Vitruvian Handle
B 8COB-426L	4 x 260mm	Vitruvian Handle

### Las Vegas

B 8LV-326-SYR	3 x 260mm	Toomey Hub
B 8LV-415L	4 x 150mm	Vitruvian Handle

### Spatula

B 8SPA-315L	3 x 150mm	Vitruvian Handle
B 8SPA-315-SYR	3 x 150mm	Toomey Hub
B 8SPA-415L	4 x 150mm	Vitruvian Handle
B 8SPA-415-SYR	4 x 150mm	Toomey Hub



### Toomey Syringe Luer Lock Adaptor - Reuseable

Convert the Toomey syringe with this reuseable, autoclavable adaptor to connect Luer Lock Fat Harvesting Cannula.

Order Ref; 1340-TOM-LLN-AL



### Toomey Syringe

Sterile 50/60ml Toomey style syringe to fit with Toomey style cannula. Available in a box of 5 syringes.

Order Ref; MSS160TP

Order Ref; 89-SYR-60TT Monoject syringe



### Syringe Snap Lock

Locks to help maintain syringe pressure, these “snap” locks connect to the back finger rest of the syringe case. Lock is maintained by twisting the plunger of the syringe to hold in place. Reuseable and autoclavable, these are available for 10ml and 60ml syringes and supplied individually.

Order Ref; 1290-TL10cc 10ml Syringe Lock

Order Ref; 1280-TL60cc 60ml Syringe Lock




### Re-useable Snapper

Order Ref; 098741 60/50 ml

Order Ref; 098742 20 ml


Order Ref; 098743 10 ml



### Syringe Racks Reuseable

Reuseable and autoclavable steel syringe racks for holding syringes of harvested fat prior to injection. Available in the following sizes:

Order Ref;	Rack	Syringe Size Holes
POR6383	3 x 6	Syringes 1 - 10ml
POR36113	3 x 6	Syringes 1 - 50ml



## Single Use Syringe Racks

Single Use steel syringe racks for holding syringes of harvested fat prior to injection.

Order Ref; SRFT2425 24 Hole Rack



## Luer Lock Syringe Cap

Fits all standard luer lock syringes.

Order Ref; MC1711 (box of 100)



## Luer Lock Transfer Adaptor

Male to male luer lock syringe connector to transfer fat from one syringe to another.


Order Ref; 993922T (box of 100)



## Hydraulic Fat Separator

Male luer lock to luer lock fat transfer adaptor, with internal reducer which changes the internal diameter of the luer from 2,5mm to 1,4mm – ideal for fat cell disruption.

Order Ref; ACC005



## Nanocube Kit

The Nanocube kit has been created as an all-in-one kit for all small to moderate volume fat grafting / transfer procedures. The Nanocube kit contains the following;

**1 x Fat harvesting cannula**

**1 x Empty IV bag** (to collect the fat)

**1 x IV tubing**

**1 x Non return 3-way tap**  
(to move the harvested fat into the IV bag in a close manor)

**1 x Nanocube**  
(use to turn MacroFat into MicroFat, Emulsified Fat or even Nanofat within the single cube – NO centrifuge required).

**1 x Micro-needle fat applicator**  
(used to place Nanofat into the epidermis)

**1 x Micro cannula**  
(for placement of microfat into the sub-dermal areas).

Order Ref; NANOKIT



### Variable Speed High Flow Infiltration Pump – with LED Fluid Volume Read-Out

Vitruvian Infiltration pump can be used in conjunction with Aquavage, Soft Touch irrigation tubing and Vibrasat for a complete closed and Power Assisted Fat grafting / body contouring procedure (as described by Del Vecchio & Wall)

- Reverse switch for flow direction preference
- Flow rate adjustable from; 0ml – 25ml - 475ml per minute
- Display: flow rate and volume dispensed
- Advanced technology allows for fluid volume measurements
- Backflow feature minimises fluid leakage
- Pneumatic foot switch activated
- Option to change from “momentary” or “maintained” foot control

Order Ref; B89055

B89213 - Soft touch tubing - for use with the B89055 pump



### K Pump Infiltration

The K Pump is designed to infiltrate tumescent fluid into a patient by peristalsis. It has a footswitch to activate it, and the speed of infiltration is adjustable

Order Ref; KPump



### Liposat® power Infiltration Device

Liposuction using the tumescent local anesthesia (TLA) technique is the most frequently performed procedure around the world. The infiltration pump Liposat® power has been specially developed for this application. Its features and functions have been designed and adjusted for surgical and aesthetic medicine. This method has principal advantages for any treatment that requires local anaesthesia, such as liposuction or vascular surgery. Liposat Infiltration pump can be used in conjunction with Aquavage, Soft Touch irrigation tubing and Vibrasat for a complete closed and Power Assisted Fat grafting / body contouring procedure (as described by Del Vecchio & Wall).

Order Ref; 00002274

## Tumescent Infiltration Tubing

### Tumescent Infiltration Tubing

High flow Irrigation pump tubing suitable for rotary motor irrigation.

Order Ref; TP01 - for narrow (European) gauge irrigation pumps.  
(Single spike, box of 10).

B89213 - Soft Touch Tubing – for infiltration of fluids and fat in closed system from Aquavage cansiter.


B89216 - 15' Large bore tubing with silicone insert for infiltration pumps (box of 10)





## Tumescence Tube for Liposat® power

Order Ref; 00002251



## Infiltration Cannulae - Single Use

Sterile, single-use infiltration/irrigation cannulae with Luer Lock connector and multiple spiral hole tips. The codes with a reference number beginning ABME all have metal luer lock hubs. Available in the following sizes:

Order Ref;	ø (mm)	Length (mm)	Tip Holes	
CI16150	1.65	150	10	Short - ideal for neck & arms
CI18100	1.25	100	10	
CI21100	2.10	100	10	
CI21150	2.10	150	10	
CI21250	2.10	250	10	
CI25150	2.50	150	10	Long - ideal for trunk
CI25250	2.50	250	10	
CI25300	2.50	300	10	
ABME25150	2.50	150	10	Metal Hub
ABME1143	2.50	200	10	
ABME1144	2.50	300	10	



## Infiltration Handles – Reusable

Infiltration handle with on/off control and lock. Luer lock fitting at both ends allow for connection to Infiltration tubing (see TP01) and Infiltration cannula (see below)

Order Ref;	B89TIH-C	With on/off control
	B89IHN-C	Without control



## Vitruvian Infiltration Cannulae - Re-usable

Order Ref;	ø (ga)	Length (mm)
Body:		
B 8TI-1015-LL	10	150
B 8TI-1032-LL	10	320
B 8TI-1215-LL	12	150
B 8TI-1220-LL	12	200
B 8TI-1226-LL	12	260
B 8TI-1232-LL	12	320
B 8TI-1236-LL	12	360
B 8TI-1415-LL	14	150
B 8TI-1420-LL	14	200
B 8TI-1426-LL	14	260
B 8TI-1432-LL	14	320
B 8TI-1436-LL	14	360
B 8TI-1615-LL	16	150
B 8TI-1620-LL	16	200
B 8TI-1626-LL	16	260
B 8TI-1632-LL	16	320
Facial:		
B 8FAN-1415-LL	14	150
B 8FAN-1615-LL	16	150





### Macro-Fat Grafting Kit (Ref MACK02)

Single-use Macro-Fat Grafting Kit of autologous adipose tissues.

Comprising of:

1 x Klein Cannula for infiltration - [Order Ref; CI25150](#)

- 2.5mmØ x 150mm with 10 exit holes and Luer lock hub

1 x Autologous Fat Harvesting Cannula - [Order Ref; PICK3518](#)

- 3.5mm Ø x 180mm, 5 side holes and luer lock hub

2 x Transfer Luer Lock Connector - [Order Ref; 993922T](#)

1 x Fat injection cannula – [Order Ref; F9020100](#)


- 2.0mm Ø x 100mm with 1 side exit, luer lock hub

1 x Fat Injection Cannula – [Order Ref; F9015100](#)

- 1.5mm Ø x 100mm with 1 side exit, luer lock hub

4 x Syringe bung /stand – [Order Ref; MC1711](#)

- 4 syringe caps that can be used to stand the harvested fat prior to separation of fluid portion



### Micro-Fat Grafting Kit (Ref MFGK02)

Single-use Fat Grafting Kit for Micro-Fat Grafting of autologous adipose tissues.

Comprising of:

1 x Hydraulic fat separator with reduction - [Order Ref; ACC005](#)

- from 2.5-1.4mm luer lock to luer lock connector

1 x Micro-Fat injection cannula – [Order Ref; F9010070](#)

- 20g (1.0mm) x 70mm single hole exit with luer lock for Micro-fat injection

1 x Micro-Fat injection cannula – [Order Ref; MI2250KR](#)

- 22g (0.71mm) x 50mm single hole exit with luer lock for Nano-fat injection

1 x Tonnard fat harverst – [Order Ref; TSFH2003](#)

- 200mm x 3mm x 24 hole


1 x Fat transfer connector – [Order Ref; 993922T](#)



### Nanocube

Single-use closed sterile system for mechanical processing of lipo aspirate (autologous fat tissue) into milli, micro and nano fat grafts according to the depth of transfer.

[Order Ref; Nanocube](#)



## AquaVage – High Volume Fat Grafting

AquaVage is a sterile, single-use, closed system for fat harvesting. AquaVage is a container, which can be connected in-line with any aspiration pump. The fat is harvested using your preferred liposuction cannula at reduced flow volumes (18 inHg/500mmHg).

The AquaVage container separates the fat and fluids which allows the surgeon to first remove the fluid portion and then the fat cells from the container via a sterile and closed method.

A sterile 60cc toomey syringe is included along with the canister hoses to facilitate easy removal of the fat and transfer into Luer Lock syringes ready for reinjection.

AquaVage offers these major advantages over other fat harvesting systems;

- Can be used with all existing equipment, including wall suction and all types of Liposuction cannula and equipment, including Power Assisted and Ultrasonic Liposuction.

- Closed sterile, single-use system
- No Centrifuge required
- Excellent quality fat, no blocking even with fibrous fat
- Easy to use and operate
- Can be used as a completely closed system for higher volume fat grafting with either the Liposat, or the Black and Black infiltration pump with the B89213 soft touch tubing.

Each AquaVage kit contains a sterile Toomey Syringe, outlet tubing with connectors for Toomey and Luer lock and a sterile suction tubing set to connect to the aspiration device.

Order Ref;

AV1200

AquaVage 1,2 litre

AV2000

AquaVage 2,0 litre

AVSCS

AquaVage Stand (reuseable)



## Luer Lock Liposuction / Fat Harvesting Cannulae

An extensive range of excellent quality sterile, single-use cannulae that are packaged individually and available in a box of 10. All cannulae are made from medical grade highly polished steel, both internal and external, to prevent or reduce blocking of the cannulae during harvesting of the fat and also to reduce friction injury to the skin.

Cannulae are available in a wide variety of widths (Ø mm) from 1.45mm diameter or 17g through to 3.50mm or 11g and in a wide variety of lengths from 85mm to 300mm. Cannulae have a wide variance of tip styles, holes and layout as described in the pictures. All cannulae are fitted with a Luer Lock standard hub that can be used with all Luer Lock Syringes and LipiVage.

		Order Ref;	Holes	Ø (mm)	Length (mm)
MACRO FAT Ideal for breast & buttocks		PICK3518	5	3.5	180
		CP32301	1	3,00	230
		CP42301		4,00	230
		CP3015	2	3,00	150
		CP3023		3,00	230
		CP17853HS	3	1,45	85
		CP161053HS		1,65	105
		CP141253HS		2,10	125
		CP30153S		3,00	150
		CP141253HV		2,10	125
		CP31503M	MERCEDES	3,00	150
		CP32303M		3,00	230
		CP42303M		4,00	230
		CP32303T		3,00	230
		CP42303T		4,00	230
		CP141304H	4	2,10	130
		CP30154H		3,00	150
		CP141256HMP	6	2,10	125
		CP141506HS		2,10	150
		CP35176H		2,50	170
		CP122007H	7	2,80	200
		CP123007H		2,80	300
		CP141508H	8	3,00	150
MICRO FAT		CP301512R	12	3,00	150
		CP301524H	24	3,00	150
		CP302324H		3,00	150









## Cannulae for Injection

Precision made cannula for injection of fluids, dermal fillers and fat. Made from Medical grade steel, polished both internally and externally, to reduce cannula blocking and damage to the skin due to friction. Cannulae cause less bleeding (haematoma) and injury to capillaries and result in less trauma and pain when used, compared to needles.

## Micro-Cannulae

Available in a range of sizes from 18g (1.24mm) to 30g (0.31mm) and a variety of cannulae length, with Luer Lock hub.

Box of 20 units, packaged sterile

	REF.	Ø (mm)	Longueur Length (mm)
	MI18110KR	18G / 1.2mm	110
	MI1890KR		90
	MI1870KR		70
	MI2170KR	21G / 0.81mm	70
	MI2140KR		40
	MI2290KR	22G / 0.70mm	90
	MI2270KR		70
	MI2250KR		50
	MI2370KR	23G / 0.63mm	70
	MI2350KR		50
	MI2450KR	24G / 0.55mm	50
	MI2438KR		38
	MI2550KR	25G / 0.50mm	50
	MI2650KR	26G / 0.45mm	50
	MI2638KR		38
	MI2625KR		25
	MI2740KR	27G / 0.40mm	40
	FILL0860	21G / 0.80mm	60
	FILL0440	27G / 0.40mm	40






## Fat Injection Cannulae – Single-Use




A range of precision made cannulae, with polished steel, both externally and internally to prevent the clogging and blocking of the cannula during injection of more viscous material, such as fat.

These cannulae are available with a Luer Lock hub and individually sterile, packaged in a box of 10.



Example of 45°cannula hub

		REF.	Hole	Ø(mm)	Length (mm)
FACIAL / NECK		CF0965290	2	0,90	65
		CFL201502LE		2,10	150
		CF15745	1 / 45°	1,65	70
		CF151045		1,65	100
		CF201545		2,10	150
		CF202345		2,10	230
		CF10790	1 / 90°	1,05	70
		CF12790		1,25	70
		CF14790		1,45	70
		CF15790		1,65	70
CF16890		1,65		80	
CF151090		1,65		100	
CF201090		2,10		100	
CF201590		2,10		150	
CF202290		2,10		220	
CF27990		2,75		90	
CF271590		2,75		150	
CF301590		3,00		150	
CF301990		3,00		190	
CF401990		4,00		190	
SCAR RELEASE CANNULA		CP1690V	1	1,65	90
		CP16150V		1,65	150
		CP1490V		2,10	90
		CP14150V		2,10	150
		F9010070	1	1,0	70
		F9012070		1,2	70
		F9015070		1,5	70
		F9015100		1,5	100
F9020100		2,0		100	
F9027150		2,7		150	


	Order Ref;	Ø (mm)	Length (mm)	Shape
 <p>Style I</p>	B8FIC-I-1415-LL	14	150	Straight
	B8FIC-I-1420-LL	14	200	Straight
	B8FIC-I-1425-LL	14	250	Straight
	B8FIC-I-167-LL	16	70	Straight
	B8FIC-I-169-LL	16	90	Straight
	B8FIC-I-187-LL	18	70	Straight
	B8FIC-I-189-LL	18	90	Straight
	B8FIC-I-195-LL	19	50	Straight
	B8FIC-I-CC1415-LL	14	150	Concave
	B8FIC-I-CC1610-LL	16	100	Concave
	B8FIC-I-CC169-LL	16	90	Concave
	B8FIC-I-CV1415-LL	14	150	Convex
	B8FIC-I-CV1610-LL	16	100	Convex
	B8FIC-I-CV169-LL	16	90	Convex
 <p>Style II</p>	B8FIC-II-167-LL	16	70	Straight
	B8FIC-II-169-LL	16	90	Straight
	B8FIC-II-187-LL	18	70	Straight
	B8FIC-II-189-LL	18	90	Straight
	B8FIC-II-CC1610-LL	16	100	Concave
	B8FIC-II-CC169-LL	16	90	Concave
	B8FIC-II-CV1610-LL	16	100	Convex
	B8FIC-II-CV169-LL	16	90	Convex
 <p>'V' Style</p>	B8FIC-V-1215-LL	12	150	Straight
	B8FIC-V-1225-LL	12	250	Straight
	B8FIC-V-1415-LL	14	150	Straight
	B8FIC-V-1425-LL	14	250	Straight

## Female Garments




### 9" Abdominal Binder

Using a unique cotton fabric, 9" (23cm) panel structure, with a full length Velcro fastening, the binder can be fitted to all sizes and the amount of compression adjusted easily and on a daily base, that will allow you to remain in comfort whilst you assist and recover from surgery. - [Order Ref; 355-9](#)




### Full body Girdle - Zippered

Stage 1 Zippered Garments allow greater ease when getting the garment on / off particularly important in the early stages of recovery after surgery. The Stage 1 garments are designed to allow for greater expansion to accommodate post-surgical swelling with appropriate compression to match. - [Order Ref; 855](#)



### Abdominal Girdle - Zippered

The Female Abdominal Girdle Compression Cosmetic Surgery Garment (with Zipper) – Stage One features zippers for easy application, adjustable shoulder straps, targeted compression for the abdominal area, reinforced flat seams, gripper lace to prevent roll up, and a reinforced open crotch. - [Order Ref; 850](#)




### Above-knee Girdle - Zippered

The Female Above Knee Girdle Compression Garment (with Zipper) – Stage One Garment improves blood circulation, minimises swelling after the procedure, flushes the body out of potentially harmful fluids, accelerates the healing process, and allows the patient to return to daily routines sooner. - [Order Ref; 753](#)



### Abdominal Brief - Zippered

The Female Abdominal Brief Cosmetic Surgery Compression Garment (with Zipper) – Stage One Garment improves blood circulation, minimizes swelling after the procedure, flushes the body out of potentially harmful fluids, accelerates the healing process, and allows the patient to return to daily routines sooner. - [Order Ref; 75](#)



### Compression Arm Sleeve With Adjustable Cotton Knit Bra


Made with 3/4 inch length compression arm sleeves, built-in seamless cotton knit bra, convenient front and back hook and eye closure, and no wiring and no seams at the surgical areas for maximum comfort and to avoid irritation, the Plastic Surgery Compression Arm Sleeves w/Adjustable Cotton Knit Bra – Moderate Support is ideal for a patient having multiple procedures. - [Order Ref; B943](#)

## Male Garments




### 12" Triple Panel Binder

This 12" triple panel binder is perfect following abdominal and flank procedures. It is made from strong, durable stretch fabric which is lightweight, cool and comfortable. - [Order Ref; 340-12](#)




### Zippered Abdominal / Chest Garment

The Zippered Abdominal/Chest Garment is ideal for chest, abdominal & flank procedures. The breathable powernet fabric retains compressive memory is washable and fast drying and is lightweight, cool and comfortable. - [Order Ref; V840](#)



### Zippered Above-knee Abdominal Garment

This garment is ideal following abdominal, flank, hip & thigh procedures. It has been carefully designed to incorporate side zip with cotton guard tape & high waist design with 3" wide waistband for ultimate comfort post operation. - [Order Ref; 743](#)



### Adjustable Compression Vest

The Male Cosmetic Surgery Compression Vest provides compression in the abdomen, chest, and flank areas and is appropriate for liposuction, gynecomastia, and pectoral implants. This garment improves blood circulation, minimizes swelling after the procedure, flushes the body out of potentially harmful fluids, accelerates the healing process, and allows the patient to return to daily routines sooner. - [Order Ref; 640](#)

## Accessories



### Bruise MD™ 30ml

Bruise MD™ is a unique formulation using ingredients that are tried and tested to aid bruise prevention and help improve existing bruises. - [Order Ref; BR2015](#)



### Epifoam Compression Pads

Epifoam offers patients an important evolution in liposuction recovery. Providing a superior platform over which to apply post-operative compression garments. - [Order Ref; EPF30](#)



### SILDERM™ Dual Action Scar Gel

SilDerm™ Dual Action Scar Gel is a lightweight, self-drying silicone gel for the treatment of scars. - [Order Ref; 387-0961](#)

## Our Garments Range



### Female Garments

#### Abdominal & Mid Body

- 9" Abdominal Binder - [Order Ref; 355-9](#)
- 6" Single Panel Binder - [Order Ref; 350-6](#)
- 9" Double Panel Binder - [Order Ref; 350-9](#)
- Torso/Brief Zippered - Stage 1 - [Order Ref; 85](#)
- Torso/Brief Non-Zip - Stage 2 - [Order Ref; 165](#)
- High back Torso/Brief With Bra Zippered - Stage 1 - [Order Ref; B85](#)
- High back Torso/Brief With Bra Non-Zip - Stage 2 - [Order Ref; B165](#)

#### Lower Body

- Labia Support Thong - [Order Ref; 1052](#)
- Labia Support Brief - [Order Ref; 1054](#)
- Below-knee Girdle Zippered - Stage 1 - [Order Ref; 755](#)
- Below-knee Girdle Non-Zip - Stage 2 - [Order Ref; 1455](#)
- Above-knee Girdle Zippered - Stage 1 - [Order Ref; 753](#)
- Above-knee Girdle Non-Zip - Stage 2 - [Order Ref; 1453](#)
- Tram Flap Girdle Zippered - Stage 1 - [Order Ref; 750](#)
- Tram Flap Girdle Non-Zip - Stage 2 - [Order Ref; 1450](#)

#### Upper Body

- Dorsocervical Female Garment - [Order Ref; DC-65](#)
- Long Arm Sleeve With Lace - [Order Ref; L933](#)
- Compression Arm Sleeve With Adjustable Bodice - [Order Ref; 943](#)
- Compression Arm Sleeve With Adjustable Bra - [Order Ref; B943](#)
- Long Arm Sleeve Bolero - [Order Ref; BV933](#)
- Short Arm Sleeve Bolero - [Order Ref; BV930](#)

#### Full Body

- Full Body Girdle Zippered - Stage 1 - [Order Ref; 855](#)
- Full Body Girdle Non-Zip - Stage 2 - [Order Ref; 1655](#)
- High-Back Body Girdle Zippered - Stage 1 - [Order Ref; 853-H](#)
- High-Back Body Girdle Non-Zip - Stage 2 - [Order Ref; 1653-H](#)
- Below-knee Body Girdle Zippered - Stage 1 - [Order Ref; BA855](#)
- Below-knee Body Girdle Non-Zip - Stage 2 - [Order Ref; B1655](#)
- Above-knee Body Girdle Zippered - Stage 1 - [Order Ref; B853](#)
- Above-knee Body Girdle Non-Zip - Stage 2 - [Order Ref; B1653](#)

### Male Garments

#### Abdominal

- 9" Abdominal Binder - [Order Ref; 345-9](#)
- 12" Triple Panel Binder - [Order Ref; 340-12](#)
- 9" Double Panel Binder - [Order Ref; 340-9](#)

#### Upper Body



- Dorsocervical Male Vest - [Order Ref; DC-64](#)
- Zippered Compression Vest With Arms - Stage 1 - [Order Ref; 642](#)
- Non-zip Compression Tank Top - Stage 2 - [Order Ref; 1241](#)
- Zippered Compression Tank Top - Stage 1 - [Order Ref; 641](#)
- Compression Vest - Stage 2 - [Order Ref; 1240](#)
- Adjustable Compression Vest - Stage 1 - [Order Ref; 640](#)

#### Lower Body

- Scrotum Support Brief - [Order Ref; 1045](#)
- Zippered Above-knee Abdominal Garment - Stage 1 - [Order Ref; 753](#)
- Zippered Gluteus Garment - Stage 1 - [Order Ref; G740](#)

#### Full Body

- Non-zip Abdominal Garment - Stage 2 - [Order Ref; 1640](#)
- Zippered Abdominal Garment - Stage 1 - [Order Ref; 840](#)
- Non-zip Abdominal/Chest Garment - Stage 2 - [Order Ref; V1640](#)
- Zippered Abdominal/Chest Garment - Stage 1 - [Order Ref; V840](#)
- Zippered Bodysuit - Stage 1 - [Order Ref; BS843](#)







Merrow Business Park,  
Guildford, Surrey GU4 7WA

T: +44 (0)1483 456 007  
E: [sales@eurosurgecal.co.uk](mailto:sales@eurosurgecal.co.uk)

[www.eurosurgecal.co.uk](http://www.eurosurgecal.co.uk)