
Stainless Steel

Tebbetts™ Platform

Grasping forceps with unique jaw configuration for maximum gripping efficiency. Facilitates harvesting of large grafts without fragmenting.

88-0677 7-1/4" 18.5 cm Angled



Stainless Steel

Sheen-Killian

Allows positive grip of elusive tissues.

88-0675 6-3/4" 17 cm Angled



Stainless Steel

Sheen Fragment

Designed for retrieving fragments. Long, thin shaft allows extension without obstructing view.

88-0085 7" 18 cm



Stainless Steel

Jansen Middleton

Large, fenestrated jaws.

88-0852 7-3/4" 19.5 cm




Fine Instruments and Video Equipment

Stainless Steel

Hartman Mosquito

Fine hemostat for delicate tissues.

88-0300	3-1/2"	9 cm	Straight
88-0301	3-1/2"	9 cm	Curved




Stainless Steel

Halstead Mosquito

Fine hemostat for plastic surgery.

88-0302	5"	12.5 cm	Straight
88-0303	5"	12.5 cm	Curved
88-0304	5"	12.5 cm	Straight 1x2 teeth
88-0305	5"	12.5 cm	Curved 1x2 teeth




Diamond-Lite®

Titanium Mosquito

Titanium for its light weight and sensitivity. The ultimate hemostatic forcep.

32-4300	3-1/2"	9 cm	Straight
32-4301	3-1/2"	9 cm	Curved
32-4302	5"	12.5 cm	Straight
32-4303	5"	12.5 cm	Curved




Stainless Steel

Delicate Mosquito

Light-weight mosquito forceps with fine jaw, ideal for grasping delicate tissue.

88-0307	4-1/2"	11.5 cm	Short jaw	Curved
88-0308	4-1/2"	11.5 cm	Regular jaw	Straight
88-0309	4-1/2"	11.5 cm	Regular jaw	Curved
88-0310	5-3/4"	14.5 cm	Regular jaw	Straight
88-0311	5-3/4"	14.5 cm	Regular jaw	Curved



Stainless Steel

Kelly

Features serrations over half the length of the jaws.

88-0313	5-1/2"	14 cm	Straight
88-0314	5-1/2"	14 cm	Curved



Stainless Steel

Crile

Popular general purpose hemostat which features serrations over the entire jaw.

88-0317	5-1/2"	14 cm	Curved
88-0319	6-1/4"	16 cm	Curved



Stainless Steel

Kocher

Heavy jaws with 1x2 teeth.

88-0330	5-1/2"	14 cm	Straight
88-0334	6-1/4"	16 cm	Straight



Stainless Steel

Adson

Long hemostat with strong curve.

88-0358	7-1/4"	18.5 cm	Curved
---------	--------	---------	--------



Fine Instruments and Video Equipment

Stainless Steel

Julian Mosquito

Extremely delicate hemostat. Sensitive to touch.

88-0359	7-1/2"	19 cm	Curved
88-0360	9-1/4"	23.5 cm	Curved



Stainless Steel

Lahey Vulsellum

Provides secure grip of undermined flaps.

88-0363	6"	15.5 cm
---------	----	---------



Stainless Steel

Allis

Delicate teeth provide secure atraumatic grip.

88-0349	5-1/2"	14 cm
88-0350	6"	15 cm



Stainless Steel

Backhaus Towel

Sharp, pointed tips overlap to ensure maximum holding ability.

88-0342	3-1/2"	9 cm
88-0343	5-1/2"	14 cm



Stainless Steel

Edna Towel

Non-perforating jaws.

88-0344 4" 10 cm

88-0345 5-1/2" 14 cm



Stainless Steel

Foerster Sponge-Holding Forceps

Highest quality stainless steel.

88-0346 7" 18 cm

88-0347 9-1/2" 24 cm



Fine Instruments and Video Equipment

Stainless Steel

Tebbetts™

Designed for precise resections. Fine tip affords maximum visibility and precision. Handle and spring designed for optimal cutting efficiency and comfort.

88-3015 6-1/4" 16 cm 2.5 mm jaw



Stainless Steel

Beyer

Delicate rongeur with double-action for increased strength and ease of cut.

88-3010 7" 18 cm 3 mm jaw



Stainless Steel

Takahashi Cutting Forceps

Long, slender shaft allows removal of small fragments or polyps.

88-3050 5-1/2" 14.5 cm Shaft length 2.5 mm bite



Stainless Steel

Kleinert-Kutz Forceps

88-3012 5-1/2" 14 cm Straight

88-3013 5-1/2" 14 cm Curved

