



# EuroSurgical

## Liposuction and Lipo-modelling



## Liposuction and Lipo-modelling



# EuroSurgical


Merrow Business Park, Guildford, Surrey GU4 7WA  
T: +44 (0)1483 456007 F: +44 (0)1483 456008  
E: [info@euroSurgical.co.uk](mailto:info@euroSurgical.co.uk)  
W: [www.euroSurgical.co.uk](http://www.euroSurgical.co.uk)



# EuroSurgical


# Single use disposable Liposuction Cannula Range

We offer a comprehensive range of disposable equipment for plastic, aesthetic and cosmetic surgeons that is safe, effective and economic. The MAC III, Mercedes II & Illouz require a handle Ref: DOT0155




## Mac III

Reference	Diameter	Length
DOMA3150	3MM	150MM
DOMA3250	3MM	250MM
DOMA4150	4MM	150MM
DOMA4250	4MM	250MM




## Mercedes II

Reference	Diameter	Length
DOME3150	3MM	150MM
DOME3250	3MM	250MM
DOME4150	4MM	150MM
DOME4250	4MM	250MM



## Illouz

Reference	Diameter	Length
DOIL3150	3MM	150MM
DOIL3250	3MM	250MM
DOIL4150	4MM	150MM
DOIL4250	4MM	250MM
DOIL5250	5MM	250MM



## Disposable Liposuction Cannula

Reference	Diameter	Length	Holes
LC 31302	3.0	130	2
LC 31303	3.0	130	3
LC 313024H	3.0	130	24
LC 31304H	3.0	130	4
LC 31503HVS	3.0	150	3
LC 32302	3.0	230	2
LC 32303	3.0	280	24
LC 32303HVS	3.0	280	24
LC 32804H	4.0	230	2
LC 328024H	4.0	230	3
LC 42302	4.0	280	2
LC 42303	4.0	280	3
LC 428024H	4.0	280	24
LC 52802	5.0	280	2
LC 52803	5.0	280	3
LC 62802	6.0	280	2
LC 62803	6.0	280	3

# Liposuction Accessories



## Syringe Stop Lock and Tommey String

1280-TL60CC Toomey 50/60cc

1290-TL10CC BD 10cc



## Medela Disposable Suction Bags

077.0084 Disposable 1.5 ltr suction bag with solidifier

077.0087 Disposable 2.5 litre suction bag with solidifier



## Infiltration Tubing

Infiltration Tubing. Double spike, to use with high power assisted pump fluid drive.

ATP 845  
Box of 20 units



## Suction Tubing

ST 30200 - 10	200CM
1201HPF 1780	400CM
- Reinforced high power pressure suction tubing	


Box of 10 units



## Disposable Infiltration Cannula

Ref	Length	(dia.mm)	Holes
DON1513	150	1.3	5
DOIN2525	250	2.5	5
CI25150	150	2.5	10
C125300	300	2.5	10

# Lipo - Modelling Lipivage



An Innovation in Fat Harvesting...

LipiVage™ is the ideal tool for harvesting the ideal filler.

Many agree that fat makes an ideal dermal filler. It is natural, readily available and may even provide permanent results. However, until now fat harvesting has required relatively difficult, inefficient and time-consuming methods.

In addition, traditional methods of processing fat may have been detrimental to cell viability and lasting results.

LipiVage™ was specifically created for harvesting fat simply, quickly and carefully giving fat a better chance for survival.

## Lipivage

**D11-0290-102 standard volume**  
30ML Syringe with internal filter, for fast effective fat harvesting single use.

## New High Volume

**D11-0290-106**  
60ML Syringe with wide bore Luerlock, Complete with fat harvesting Cannula.



# Lipo - Modelling Aquavage



The unique simple design of the AquaVage gives you these advantages:

- Fast and Trouble Free
- Closed, Disposable Sterile System
- Separates Fat Automatically
- Available in 1 or 2 litre canisters
- No clogging, Even in High Fibrous Tissue
- No Centrifuge Required
- Works with any controlled aspirator
- True High Volume System

## Aquavage

AV1200 1200CC  
AV2000 2000CC  
Reusable stand AVSCS

## Single use, Disposable Fat Harvest and Fat Injection Cannula

### Micro Cannula for Filler injection

Length	110mm	90mm	70mm	50mm	40mm	38mm	25mm
18G 1,24mm	M118110		M11870				
21G 0,81mm			M12170KR				
22G 0,71mm		M12290KR	M12270KR	M12250KR			
23G 0,63mm			M12370KR	M12350KR			
24G 0,56mm				M12450KR		M124338KR	
25G 0,51mm				M12550			
26G 0,46mm				M12650KR		M12638KR	M126245KR
27G 0,41mm					M12740		
30G 0,31MM							M133025

Packaging: Box of 20, sterile



### Fat Harvest Cannula

Reference	Diameter	Length	Holes
CP 30154H	3.00	150mm	4
CP 301524H	3.00	230mm	24
CP3015	3.00	150mm	2
CP3023	3.00	230mm	2



### Fat Injection Cannula

Reference	Diameter	Length
F9010070	1,0MM	70MM
F9012070	1,2MM	70MM
F9015100	1,5MM	100MM
F9015070	1,5MM	70MM
F9020100	2,0MM	100MM
F9027150	2,7MM	150MM



### Fat Harvest Cannula

Reference	Diameter	Length	Holes
PICU3518	3.5MM	180MM	5

## Disposable Lipofilling Kit



### Standard Kit

Order code: KLF004 KIT 4

Sterile single-use, disposable self contained kit for fat harvesting and transfer. The kit contains the following:

- 4 x 10ml Syringes, Luer lock
- 4 x Luer lock caps
- 3 x 1ml syringe, luer lock
- 1 x Adaptor, luer lock to luer lock
- 1 x Cannula 16g x 80mm, fat harvesting
- 1 x Cannula 14g x 125mm, fat re-injection
- 1 x Needle 16g x 25mm
- 1 x Needle 21g x 40mm
- 1 x Needle 30g x 13mm

### Advanced Kit

Order code: KLF003 KIT 3

Sterile single-use, disposable self contained kit for fat harvesting and transfer. The kit contains the following:

- 4 x 10ml Syringes, Luer lock
- 4 x Luer lock caps
- 4 x 1ml syringe, luer lock
- 1 x Adaptor, luer lock to luer lock
- 1 x Infiltration cannula 3 x 150mm
- 1 x Coleman type fat harvesting cannula 1,5 x 150mm
- 1 x Cannula 17g x 80mm, fat harvesting
- 1 x Cannula 18g x 70mm, fat re-injection
- 1 x Needle 16g x 25mm
- 1 x Needle 21g x 40mm
- 1 x Needle 30g x 13mm