

Footswitches

For Electrosurgery Generators







Easy
Clean

This new range of Prima "Easyclean" electrosurgery footswitches is produced in a special latex free, synthetic rubber - easy to clean!

Non-slip and water tight, the "Easyclean" range of footswitches are designed to resist the most caustic cleaning agents used in hospitals today. Your guarantee of safety and reliability.

- **Easy to clean**
- **Non slip**
- **Latex free**
- **Resists caustic cleaning agents**



GENERATOR CONNECTOR		Yellow (cut)	Blue (coag)	Black	White	Twin
	 3pin Conmed			FS212	FS213	
	 4pin Conmed/Aspen					FS250
	 Bard/Davol					FS041
	 4pin Valleylab, Megadyne					FS206
	 3pin Valleylab, Megadyne Bipolar			FS210	FS211	
	 Birtcher			FS555		
	 Erbe Bipolar				FS501	
	 Erbe					FS720
	 Eschmann 800 Series	FS209	FS2117			FS370
	 Eschmann 300 & 400 Series	FS109	FS117	FS718	FS719	FS270
	 (without connectors)	Yellow (cut) Single Footswitch FS355C	Blue (coag) Single Footswitch FS355D	Black Single Footswitch FS355E	White Single Footswitch FS176	Twin Footswitch FS176

Prima Safety features



Individually labelled to identify product



Unique serial number for tracing

Easy Clean



Eurosurgical Ltd

Merrow Business Centre,
Guildford, Surrey, GU4 7WA

+44 1483 456007

sales@eurosurgical.co.uk

www.eurosurgical.co.uk

