

Nonstick Electrodes



PTFE coated electrodes

- **Frustration free** – lets you focus on the surgery
- **Focussed energy** – puts the power where you need it
- **No scratch pad required**

Control without compromise

Electrosurgery electrodes - where technology matters...

The build-up of carbonised material (eschar) on the electrosurgery tip during surgery is radically reduced with the Prima nonstick PTFE coated electrodes.

The eschar layer prohibits the electrosurgery current being passed through the electrode, necessitating higher output settings which can cause damage to the peripheral tissue.

Prima nonstick electrodes are designed to focus the energy on the target tissue, making surgical time more productive.

- Focussed electrosurgery
- Quicker procedures
- Easy cleaning of tip
- Lower power settings
- Reduction of peripheral damage

High Temperature Shielded Tip Electrodes

- New modified high temperature resistant PTFE insulation withstands temperatures up to 350 °C
- Shielded tip electrodes offer a greater level of safety when working in cavities and areas of restricted vision, preventing unwanted contact with unintended tissue.
- Improved safety

Blades



PMS 2950
Nonstick blade electrode
64mm, shielded tip
Single use, Pack of 10



PMS 2920
Nonstick blade electrode
64mm
Single use, Pack of 10



PMS 2951
Nonstick blade electrode
70mm, shielded tip
Single use, Pack of 10



PMS 2921
Nonstick blade electrode
70mm
Single use, Pack of 10



PMS 2926
Nonstick blade electrode
100mm, shielded tip
Single use, Pack of 10



PMS 2925
Nonstick blade electrode
100mm
Single use, Pack of 10



PMS 2928
Nonstick blade electrode
130mm, shielded tip
Single use, Pack of 10



PMS 2927
Nonstick blade electrode
130 mm
Single use, Pack of 10



PMS 2933
Nonstick blade electrode
165mm, shielded tip
Single use, Pack of 10



PMS 2929
Nonstick blade electrode
165mm
Single use, Pack of 10



PMS 2931
Nonstick blade electrode
200mm
Single use, Pack of 10

Nonstick Electrodes



An extra long 130mm version of this ball electrode is available by adding "L" to the end of the code.



An extra long 130mm version of this ball electrode is available by adding "L" to the end of the code.



An extra long 130mm version of this ball electrode is available by adding "L" to the end of the code.

Needles

PMS 2934

Nonstick needle electrode
152mm
Single use, Pack of 10

PMS 2935

Nonstick needle electrode
152mm, shielded tip
Single use, Pack of 10

PMS 2945

Nonstick needle electrode
100mm
Single use, Pack of 10

PMS 2946

Nonstick needle electrode
100mm, shielded tip
Single use, Pack of 10

PMS 2922

Nonstick needle electrode
70mm
Single use, Pack of 10

PMS 2972

Nonstick needle electrode
70mm, shielded tip
Single use, Pack of 10

Blunt Needles



PMS 2910

Nonstick blunt needle electrode
70mm, shielded tip
Single use, Pack of 10



PMS 2911

Nonstick blunt needle electrode
70mm
Single use, Pack of 10

Balls

PMS 2905

Nonstick ball electrode
4mm ball, 50mm long
Single use, Pack of 10

PMS 2909

Nonstick ball electrode
5mm ball, 70mm long
Single use, Pack of 10

PMS 2908

Nonstick ball electrode
3mm ball, 110mm long
Single use, Pack of 10

PMS 2907

Nonstick ball electrode
4mm ball, 110mm long
Single use, Pack of 10

PMS 2906

Nonstick ball electrode
5mm ball, 110mm long
Single use, Pack of 10

Nonstick Electrodes



Fingerswitch



PMS252 - Box of 25
Single use fingerswitch pencil with
70mm nonstick blade electrode



PMS254 - Box of 25
Single use fingerswitch pencil with
70mm shielded nonstick blade electrode



Control without compromise

Eurosurgical Ltd

Merrow Business Centre,
Guildford, Surrey, GU4 7WA

+44 1483 456007

sales@eurosurgical.co.uk

www.eurosurgical.co.uk

