



# Facial Surgery

# PELNAC Artificial Dermis



Artificial dermis PELNAC consists of two layers; a porcine tendon-derived atelocollagen sponge layer and silicone film. It is suitable for use in full-thickness skin defect wounds and used as

alternative materials for traditional treatment for the formation of new dermis-like tissue by invasion of fibroblasts into the atelocollagen sponge matrix.

## Clinical Advantage

- Provides a high survival rate of secondary skin grafts and satisfactory aesthetic results.
- Thin split-thickness skin graft is achievable and it

reduces the damage of donor site, minimizes skin sacrifice.

- Minimal contraction or pigmentation after treatment

## Product Characteristic

- Made of atelocollagen derived from porcine tendon and silicone
- The soft collagen sponge structure ensures excellent contact with the irregular wound surface.

- Various types appropriate for wound condition are available.

- Easy transportation and storage by re-freeze dried condition


## Mechanism of Action

- Fibroblasts and capillaries infiltrate into atelocollagen from the recipient matrix and surrounding tissues and form good dermis-like tissue.

- Regeneration of dermis-like tissue is clearly different in features including collagen arrangement from scar tissue.

Size	Dimensions (mm)	Package (1 box)	Fenestrated (3mm)	Fortified (3mm)	Single Layer (3mm)
3S	40 × 30	1	PN-D40030	PN-F40030	PN-S40030
SS	40 × 60	1	PN-D40060	PN-F40060	PN-S40060
S	82 × 60	1	PN-D82060	PN-F82060	PN-S82060
M	82 × 90	1	PN-D82090	PN-F82090	PN-S82090
L	82 × 120	1	PN-D82120	PN-F82120	PN-S82120
LL	120 × 240	1	PN-D120240	PN-F120240	PN-S120240
3L	200 × 240	1	PN-D200240	PN-F200240	PN-S200240

## PELNAC Artificial Dermis (Size Variations)





# Implantech Facial Implants

## Conform Extended Anatomical Chin

Evolution of all of the technology and design of our previous chin implants. Incorporates the patented Conform grid pattern to allow for easier implant insertion and to closely conform to underlying bony structure.


Cat. No	Size	A	B	C	D	E	F
CEAC-S	Small	5.0	1.0	0.43	8.0	3.5	0.2
CEAC-M	Medium	5.6	1.2	0.50	10.0	4.2	0.3
CEAC-L	Large	6.2	1.4	0.60	10.0	4.5	0.3
CEAC-XL	X-Large	6.8	1.5	0.75	11.0	5.0	0.3
CEAC-S-SIZ-NS	Sizers	Complete Set		Dimensions in centimeters			



## Extended Anatomical Chin

Full lateral arms blend into mandible. Can be used with Glasgold Wafer for additional anterior augmentation.


Cat. No	Size	A	B	C	D	E	F
EAC-S	Small	5.5	1.0	0.45	9.0	3.7	0.2
EAC-M	Medium	5.6	1.2	0.50	9.5	4.2	0.4
EAC-L	Large	6.8	1.4	0.60	11.0	4.8	0.4
EAC-XL	X-Large	6.0	1.4	0.75	11.0	5.1	0.5
EAC-S-SIZ-NS	Sizers	Complete Set		Dimensions in centimeters			



## Terino Extended Anatomical Chin

Lateral arms are thinner and broader than Extended Anatomical Chin for less palpable edges. Profile is slightly squarer for a stronger jawline. Can be used with Glasgold Wafer for additional anterior augmentation.

Cat. No	Size	A	B	C	D	E	F
TEAC-S	Small	5.7	1.0	0.38	11.0	4.8	0.3
TEAC-M	Medium	5.7	1.2	0.46	11.3	5.0	0.3
TEAC-L	Large	5.7	1.3	0.58	11.5	5.2	0.3
TEAC-XL	X-Large	5.7	1.4	0.64	11.5	5.2	0.4
TEAC-S-SIZ-NS	Sizers	Complete Set		Dimensions in centimeters			






## Terino Square Chin - Style I

Provides anterior chin augmentation with a square profile .


Cat. No	Size	A	B	C	D	E	F	G	H
TSCI-S	Small	6.6	1.3	0.60	9.8	4.1	0.4	4.0	0.70
TSCI-M	Medium	6.5	1.4	0.70	10.9	4.8	0.3	4.2	0.80
TSCI-L	Large	6.6	1.6	0.80	12.0	5.4	0.3	4.6	0.93
TSCI-XL	X-Large	6.9	1.6	1.10	12.8	5.8	0.4	4.2	1.20
TSCI-SIZ-NS	Sizers	Complete Set		Dimensions in centimeters					



## Anatomical Chin

Central chin augmentation. Short lateral extensions blend into mandible.


Cat. No	Size	A	B	C	D	E
AC-S	Small	4.7	1.2	0.60	6.0	2.1
AC-M	Medium	5.3	1.2	0.90	6.6	2.5
AC-L	Large	5.1	1.2	1.00	6.2	2.5
AC-XL	X-Large	5.3	1.2	1.20	6.6	2.8
AC-SIZ-NS	Sizers	Complete Set		Dimensions in centimeters		



## Conform Binder Submalar

Softer and more flexible than the standard silicone implant, the patented Conform grid pattern readily adapts to the underlying facial structure. Augments the submalar region for a more youthful appearance. Sold in pairs.


Cat. No	Size	A	B	C
CBSM-S	Small	4.4	1.5	0.4
CBSM-M	Medium	4.9	1.8	0.5
CBSM-L	Large	5.4	1.9	0.6
CBSM-XL	X-Large	5.9	2.1	0.7
CBSM-SIZ-NS	Sizers	Complete Set		Dimensions in centimeters



## ePTFE-Coated Extended Anatomical Chin


Full lateral arms blend into mandible .

Cat. No	Size	A	B	C	D	E	F
ePTFE-EAC-S	Small	5.5	1.0	0.45	9.0	3.7	0.2
ePTFE-EAC-M	Medium	5.6	1.2	0.50	9.5	4.2	0.4
ePTFE-EAC-L	Large	6.8	1.4	0.60	11.0	4.8	0.4
ePTFE-EAC-XL	X-Large	6.0	1.4	0.75	11.0	5.1	0.5
ePTFE-EAC-SIZ-NS	Sizers	Complete Set		Dimensions in centimeters			



# Neoveil – Mucosal Template

Usable all over the human body



Hepatobiliary & Pancreatic surgery

ESD (Endoscopic submucosal dissection)

ENT, Oral & Maxillofacial surgery

Thoracic surgery


## What's NEOVEIL™?

- 100% Polyglycolic acid
- BIOABSORBABLE
- SOFT TISSUES/ORGANS REINFORCEMENT MATERIAL

## Retention of Strength

### NEOVEIL™ degrades by hydrolysis

- Completely absorbed after about 15 weeks




## Indications

1. Reinforcement for regions which do not take a long time to heal.
2. Reinforcement for regions of soft tissues/organs in thoracic, abdominal and oral surgery.
3. Prevention of air leakage from the lung, bronchi and diaphragm, where such prevention is difficult to perform by only using sutures or surgical staplers.


**NEOVEIL™ SHEET**  
Product video



## (Size Variations)



## Advantages



**PREVENTING  
AIR LEAKAGE**



**PREVENTING  
FLUID LEAKAGE**



**PREVENTING POSTOPERATIVE  
PANCREATIC FISTULA**

### Code

NV-L-05G

NV-L-03G

NV-L-015G

NV-M-05G

NV-M-03G

NV-M-015G

NV-S-05G

NV-S-03G

### Size

100 x 100 x 0.5 mm

100 x 100 x 0.3 mm

100 x 100 x 0.15 mm

100 x 50 x 0.5 mm

100 x 50 x 0.3 mm

100 x 50 x 0.15 mm

50 x 50 x 0.5 mm


50 x 50 x 0.3 mm



1 sheet / box

# Fat Harvesting Cannula

Cannulas specifically designed to perform an accurate fat sizing during harvesting thanks to our patented diamond-shaped ports with 360°cutting edges. Macrofat Picking cannula (**REF PICK352-L**) has a luer lock hub and a handle to convert to liposuction via aspiration pump.



- Small volume
- LL syringe

- Med/Large volume
- LL syringe / Pump



Metallic  
Luer/Lock



Handle

REF

Ø mm

G

↔ mm

Holes

Box

## INFILTRATION

	<b>DOIN1510-2</b>	1.5	17	100	5	10
	<b>DOIN2015-2</b>	2	14	150	5	10
	<b>CI16150</b>	1.65	16	150	10	10
	<b>CI18100</b>	1.25	18	100	10	10
	<b>CI21100</b>	2.1	14	100	10	10



















## HARVESTING

	<b>PICK3518-2</b>	3.5	10	180	5	10
	<b>PICK352E-2</b>	2.0	14	150	8	10
	<b>PICK352L-2*</b>	3.5	10	230	18	10
	<b>PICK352A-2</b>	4.0	8	300	8	10
	<b>CP3015</b>	3	11	150	2	10
	<b>CP3023</b>	4	8	230	2	10
	<b>CP161053HS</b>	1.65	16	105	3	10
	<b>CP141253HS</b>	2.1	14	125	3	10
	<b>CP30153S</b>	3	11	150	3	10
	<b>CP31503M</b>	3	11	150	3	10
	<b>CP32303M</b>	3	11	230	3	10
	<b>CP42303M</b>	4	8	230	3	10
	<b>CP32303T</b>	3	11	230	3	10
	<b>CP42303T</b>	4	8	230	3	10
	<b>CP30154H</b>	3	11	230	4	10
	<b>CP141256HMP</b>	2.1	14	125	6	10
	<b>CP141506HS</b>	2.1	14	150	6	10
	<b>CP141508H</b>	3	11	150	8	10
	<b>CP301512R</b>	3	11	150	12	10
	<b>CP301524H</b>	3	11	150	24	10
	<b>CP302324H</b>	3	11	150	24	10
	<b>CP122007H</b>	2.8	12	200	7	10
	<b>CP123007H</b>	2.8	12	300	7	10

\*Comes with a handle which also has a luer lock connection



# Fat Injection

	REF	Ø mm	G	↔ mm	Holes	Box
FAT INJECTION						
	FIL1270S-2	3.5	10	250	1	10
	FILU2515-2	2.5	12	200	1	10
	FIL2590S-2	2.7	12	200	1	10
	F9027150-2	2.7	12	150	1	10
	FIL25C5S-2	2.0	14	150	1	10
	FIL1670S-2	2.0	14	150	1	10
	F9020100-2	2.0	14	100	1	10
	F9015100-2	1.5	17	100	1	10
	F9015070-2	1.5	17	70	1	10
	F9012100-2	1.2	18	100	1	10
	F9012070-2	1.2	18	70	1	10
	CP1690V	1.65	16	90	1	10
	CP16150V	1.65	16	150	1	10
	CP1490V	2.10	14	90	1	10
	CP14150V	2.10	14	150	1	10
	MI18110KR	1.2	18	110	1	20
	MI1890KR	1.2	18	90	1	20
	MI1870KR	1.2	18	70	1	20
	MI2170KR	0.81	21	70	1	20
	MI2140KR	0.81	21	40	1	20
	MI2290KR	0.70	22	90	1	20
	MI2270KR	0.70	22	70	1	20
	MI2250KR	0.70	22	50	1	20
	MI2370KR	0.63	23	70	1	20
	MI2350KR	0.63	23	50	1	20
	MI2450KR	0.55	24	50	1	20
	MI2438KR	0.55	24	38	1	20
	MI2550KR	0.50	25	50	1	20
	MI2650KR	0.45	26	50	1	20
	MI2638KR	0.45	26	38	1	20
	MI2625KR	0.45	26	25	1	20
	MI2740KR	0.40	27	40	1	20
	POR6383	Stainless Steel Syringe Racks			18 Hole Rack	
	993922T	Luer Lock Transfer Adaptor			Box of 100	
	MC1711A	Luer Lock Syringe Cap			Box of 100	


# st'rim®

**MICROFAT** grafting kit - **REF PICS2130-2**



SIZING



360° SUCTION



MEDIUM FLOW



	TIP	A	B1	B2	B3	C1	C2	REF	Box
st'rim	Ø ←mm	2.0 mm 130	19G 50	21G 60	21G 40	14G 40	21G 40	PICS2130-2	5

## Lipocube




**Nanocube** "All in one" Fat Grafting - 4 ports for small volume fat grafting procedures;

Sterile, single-use device for re-sizing harvested fat particles into mini, micro and nanofat grafts for specific indications;

**Supplied sterile, pack size each (1) REF L-NANO**

### INDICATIONS




Milifat (port 1-2) 2.4mm<

Deep fat compartments, Fat-Lifting, Bone onlay graft



Microfat (port 3) 1.2mm<

Superficial fat compartments



Nanofat (port 4) 600 microns <

Subdermal, Intradermal, Masotherapy, Micro-needling, topical

# Retract Micro Dissection Needle

Retract is a small scale fingerswitch that contains an ultra sharp tungsten needle. The considered ergonomics of the device offers a greater level of accuracy in microsurgical procedures, especially when working around the eyes.

The soft touch buttons are placed close to the operative site to offer optimal control, more like writing with a pen.

Retract is furthermore a safety device - the needle lives inside the handle and is only extended to use the device. It can then be fully retracted when not in use and for safe disposal.

## The Surgeon's 'ballpen'


Control is absolute, so close to the operating tip that it becomes an extension of the finger, no bending or lateral movements, perfect balance between the user and microdissection surgery.

## The Nurse's friend

Needlestick injuries are minimised with the retracting needle: press to release and twist to retract back into the safety of the handle

## Three products in one

Two needles and handle all in one product




Product REF	NPC Code		QTY
PMS8001	FDJ3629	Retract electrosurgery fingerswitch with retractable micro dissection needle	10


# Supervac - Telescopic Smoke Evacuation Fingerswitch

Supervac has been specifically designed to meet the demands of surgeons. It is slimline, lightweight and offers optimal visualisation, even when working in cavities and areas of restricted vision.


- Buttons float up and down the barrel offering telescopic function for deep cavity surgery using only one product.
- Length easily adjusts in seconds.
- 360 degree swivel on back of product and non kink tubing prevents hand fatigue
- Maximum suction output over 85 Litres per minute.



**SUPER  
VAC**



Standard 3 pin  
monopolar connector.



10/22mm smoke  
evacuation connector

Product REF	NPC Code		QTY
<b>PMS6300N</b>	<b>FDJ2388</b>	SuperVac smoke evacuation fingerswitch with non-stick blade electrode and nose cone	<b>20</b>
<b>PMS6305N</b>	<b>FDJ3658</b>	SuperVac smoke evacuation fingerswitch with modified non stick blade electrode	<b>20</b>
<b>PMS6301N</b>		SuperVac smoke evacuation fingerswitch with non-stick blade electrode and microdissection needle and nose cone	<b>20</b>

\*All products with 4 metre cable, 3 metre tubing and terminated with three pin electrosurgery plug.

## Options available with Supervac



**PMS6300N** - with standard non stick blade electrode and nose cone (Compatible with all other >70mm electrodes)

**PMS6301N** - with tungsten microdissection needle and nose cone.

**PMS6305N** - with non stick electrode, shielded 5mm tip and nose cone.

## Other compatible electrodes supplied separately ...see page 17 for full range



**Non stick needle electrodes**

**Non stick ball electrodes**

**Lletz Loops**



# Cotswold - Tungsten Micro dissection needles

## Cotswold microdissection needles

Cotswold microdissection needles offer improved cutting and coagulation performance, whilst minimising thermal spread, offering the surgeon the precision required for exceptionally fine dissection in a wide variety of procedures.

## Precise electrosurgery

The Cotswold needle tungsten tip remains sharp throughout the procedure. With a melting point greater than 3420°C, the needle will not become blunt, especially at the low power settings.

## The effective alternative

The fine tungsten tip of the Cotswold needle allows for greater accuracy to be achieved at very low power settings, in both cut and coag modes. This accuracy cannot be achieved with traditional stainless steel needles. Due to the extremely small surface area, a focused spark plasma is created at very low power settings, almost eliminating tissue drag and thermal damage, which could cause scarring.


*Single use product\**



**Safe** - Unique fingertip allows for easy insertion and extraction, without limiting vision and prevents needlestick injuries




**PMS1650 / NPC FDJ3637** - Cotswold Needle 30mm shaft, 10mm straight needle, pack of 10



**PMS1651 / NPC FDJ2405** - Cotswold Needle 30mm shaft, 20mm straight needle, pack of 10




**PMS1614 / NPC FDJ3603** - Cotswold Needle 40mm shaft, 20mm straight needle, pack of 10



**PMS1654 / NPC FDJ3638** - Cotswold Needle 55mm shaft, 10mm straight needle, pack of 10



**PMS1602 / NPC FDJ3601** - Cotswold Needle 65mm shaft, 10mm straight needle, pack of 10



**PMS1656 / NPC FDJ3609** - Cotswold Needle 75mm shaft, 10mm straight needle, pack of 10



**PMS1613 / NPC FDJ3642** - Cotswold Needle 100mm shaft, 10mm straight needle, pack of 10



**PMS1612 / NPC FDJ3652** - Cotswold Needle 130mm shaft, 10mm straight needle, pack of 10




**PMS1611 / NPC FDJ3647** - Cotswold Needle 150mm shaft, 10mm straight needle, pack of 10



**PMS1610 / NPC FDJ3602** - Cotswold Needle 180mm shaft, 10mm straight needle, pack of 10





#### 45° Angled Needles 5mm Tip

**PMS1692 / NPC FDJ3615** - Cotswold Needle 30mm shaft, 5mm 45o angled needle, pack of 10



**PMS1693 / NPC FDJ3606** - Cotswold Needle 40mm shaft, 5mm 45o angled needle, pack of 10



**PMS1694 / NPC FDJ2335** - Cotswold Needle 55mm shaft, 5mm 45o angled needle, pack of 10




**PMS1695 / NPC FDJ3610** - Cotswold Needle 75mm shaft, 5mm 45o angled needle, pack of 10



**PMS1696 / NPC FDJ3630** - Cotswold Needle 100mm shaft, 5mm 45o angled needle, pack of 10



**PMS1697 / NPC FDJ3616** - Cotswold Needle 130mm shaft, 5mm 45o angled needle, pack of 10



#### 45° Angled Needles 10mm tip

**PMS1653 / NPC FDJ3598** - Cotswold Needle 30mm shaft, 10mm 45o angled needle, pack of 10



**PMS1631 / NPC FDJ3639** - Cotswold Needle 40mm shaft, 10mm 45o angled needle, pack of 10



**PMS1655 / NPC FDJ3666** - Cotswold Needle 55mm shaft, 10mm 45o angled needle, pack of 10




**PMS1603 / NPC FDJ2387** - Cotswold Needle 65mm shaft, 10mm 45o angled needle, pack of 10



**PMS1605 / NPC FDJ3649** - Cotswold Needle 75mm shaft, 10mm 45o angled needle, pack of 10




**PMS1623 / NPC FDJ3599** - Cotswold Needle 100mm shaft, 10mm 45o angled needle, pack of 10



#### 90° Angled Needles

**PMS1660 / NPC FDJ3656** - Cotswold Needle 30mm shaft, 10mm 90o angled needle, pack of 10



**PMS1661 / NPC FDJ3669** - Cotswold Needle 40mm shaft, 10mm 90o angled needle, pack of 10



**PMS1662 / NPC FDJ3605** - Cotswold Needle 55mm shaft, 10mm 90o angled needle, pack of 10



**SUPERVAC SMOKE EVACUATION**  
Product Video



**IMPLANTECH - CHIN IMPLANT**  
Product Video



**PELNAC**  
Product Video



**LIPOCUBE NANOCUBE**  
Product Video





**Eurosurgical Limited**  
1F Merrow Business Park, Guildford,  
Surrey GU4 7WA, United Kingdom

+44 (0) 1483 456 007  
[sales@eurosurgical.co.uk](mailto:sales@eurosurgical.co.uk)  
[www.eurosurgical.co.uk](http://www.eurosurgical.co.uk)