



## Body Contouring & Fat Grafting

## Infiltration

### Variable Speed High Flow Infiltration Pump - with LED Fluid Volume Read-Out

Vitruvian Infiltration pump can be used in conjunction with Aquavage, Soft Touch irrigation tubing and Vibrasat for a complete closed and Power Assisted Fat grafting / body contouring procedure (as described by Del Vecchio & Wall)



- Reverse switch for flow direction preference
- Flow rate adjustable from; 0ml - 25ml - 475ml per minute
- Display: flow rate and volume dispensed
- Advanced technology allows for fluid volume measurements
- Backflow feature minimises fluid leakage
- Pneumatic foot switch activated
- Option to change from "momentary" or "maintained" foot control

Order Ref; B89055

B89213 - Soft touch tubing - for use with the B89055 pump

### Tumescent Infiltration Tubing

High flow Irrigation pump tubing suitable for rotary motor irrigation.

Order Ref;

TP01 - for narrow (European) gauge irrigation pumps.

(Single spike, box of 10).

B89213 - Soft Touch Tubing - for infiltration of fluids and fat in closed system



B89216 - 15' Large bore tubing with silicone insert for infiltration pumps

(box of 10)




## Infiltration Cannula

Multi-spiraled cannulas for tumescent fluid infiltration dissection prior to fat harvesting



Vortex Effect



360° Distribution

## Infiltration


REF	Ø mm	G	↔ mm	Holes	Box
<b>INFILTRATION</b>					
<b>DOIN1510-2</b>	1.5	17	100	5	10
<b>DOIN2015-2</b>	2	14	150	5	10
<b>DOIN2525-2</b>	2.5	13	250	5	10
<b>CI16150</b>	1.65	16	150	10	10
<b>CI18100</b>	1.25	18	100	10	10
<b>CI21100</b>	2.1	13	100	10	10
<b>CI21150</b>	2.1	14	150	10	10
<b>CI21250</b>	2.1	13	250	10	10
<b>CI25300</b>	2.5	13	300	10	10
<b>CI25150</b>	2.5	13	150	10	10



### Infiltration Handles - Reusable

Infiltration handle with on/off control and lock. Luer lock fitting at both ends allow for connection to Infiltration tubing and Infiltration cannula

Order Ref; **B89TIH-C** With on/off control



### Vitruvian Infiltration Cannulae - Re-usable

Order Ref;                      ø (ga)    Length (mm)

Body:

<b>B 8TI-1015-LL</b>	10	150
<b>B 8TI-1032-LL</b>	10	320
<b>B 8TI-1215-LL</b>	12	150
<b>B 8TI-1220-LL</b>	12	200
<b>B 8TI-1226-LL</b>	12	260
<b>B 8TI-1232-LL</b>	12	320
<b>B 8TI-1236-LL</b>	12	360
<b>B 8TI-1415-LL</b>	14	150
<b>B 8TI-1420-LL</b>	14	200
<b>B 8TI-1426-LL</b>	14	260
<b>B 8TI-1432-LL</b>	14	320
<b>B 8TI-1436-LL</b>	14	360
<b>B 8TI-1615-LL</b>	16	150
<b>B 8TI-1620-LL</b>	16	200
<b>B 8TI-1626-LL</b>	16	260
<b>B 8TI-1632-LL</b>	16	320

Facial:

<b>B 8FAN-1415-LL</b>	14	150
<b>B 8FAN-1615-LL</b>	16	150

## Aspiration



### Medela - Dominant Flex Suction Pump

One-touch adaptable flow rates allow for a wider variety of procedures to be carried out with just one pump. Designed for use in hospitals, clinics and medical practices, the Dominant Flex is the optimal medical suction pump for all suctioning needs. Fitted with 2 x 2,5 litre bottles and a changeover valve. Optional on/off foot switch is available.

#### Order Ref;

071.003T Dominant Flex Suction Pump

077.0571 Bacteria filter

077.0563 The Medela changeover valve with vacuum tubing

077.0085 2.5l Suction jar for single use liners

077.0913 Silicone tubing with connector

077.0521 Clamp holder



### Medela - Vario 18 Suction Pump

A multipurpose medical suction pump for hospitals and clinics that is extremely quiet, offers reliability, and is lightweight, portable and robust.

High vacuum up to -75 kPa / - 563 mmHG

Flow of 181 / min

Weight: 4.2 kg

Order Ref; VARIO 18 CD



### Medela Disposable Suction Liners

For use with ALL versions of Medela Aspiration pumps, such as Vario, and Dominant flex.

#### Order Ref;

077.0083 1.5 litre suction liner (box of 40)

077.0084 1.5 litre suction liner with solidifier (box of 30)

077.0086 2.5 litre suction liner (box of 40)

077.0087 2.5 litre suction liner with solidifier (box of 30)

### Suction Tubing

Easy to fit with two flexible connectors that will accommodate all types of push-fit cannulae and connections. Available in clear high strength plastic that prevents collapsing even at high aspiration pressure. Packaged 10 per box.

#### Order Ref;

H91390E - 10 - 200cm tubing -10mm (external) 6.2mm (internal) with 11.2mm connector (internal)

H91400E - 10 - 300cm tubing -10mm (external) 6.2mm (internal) with 11.2mm connector (internal)

1201HPF1780 - 400cm reinforced corrugated tubing - 9mm (external) 5.7mm (internal) with 10mm connector (internal)

1201HPF4817 - Lipomatic air tubing




## Liposuction Cannula










Body reshaping cannulas with sharp edges ports "guillotine" design to cut adipose tissue



<p><b>DOLFIN</b> 3 Ports</p> 	<p>Medium SUCTION EFFICIENCY (avg.)</p>		<p>260° SUCTION PRECISION</p>	
<p><b>DOLFIN PLUS</b> 7 Ports</p> 	<p>High SUCTION EFFICIENCY (avg.)</p>		<p>260° SUCTION PRECISION</p>	
<p><b>DOLFIN PRO</b> 8 Ports</p> 	<p>Very High SUCTION EFFICIENCY (avg.)</p>		<p>360° FULL SUCTION</p>	



## Single Use Cannulae

	REF	Ø mm	G	↔ mm	Holes	Box
<b>LIPOSUCTION</b>						
	<b>LC31302</b>	3	11	130	2	6
	<b>LC32302</b>	3	11	230	2	6
	<b>LC42302</b>	4	8	230	2	6
	<b>LC42802</b>	4	8	280	2	6
	<b>LC31303</b>	3	11	130	3	6
	<b>LC32303</b>	3	11	230	3	6
	<b>LC42303</b>	4	8	230	3	6
	<b>LC42803</b>	4	8	280	3	6
	<b>LC52803</b>	5	6	280	3	6
	<b>LC62803</b>	6	5	280	3	6
	<b>LC32303HVS</b>	3	11	230	3	6
	<b>LC31304H</b>	3	11	130	4	6
	<b>LC32804H</b>	3	8	280	4	6
	<b>LC313024H</b>	3	11	130	24	6
	<b>LC328024H</b>	3	11	280	24	6
	<b>LC428024H</b>	4	8	280	24	6
	<b>DOLF3150-2</b>	3	11	150	3	10
	<b>DOLF3250-2</b>	3	11	250	3	10
	<b>DOLF4150-2</b>	4	8	150	3	10
	<b>DOLF4250-2</b>	4	8	250	3	10
	<b>DOLF4350-2</b>	4	8	350	3	10
	<b>DOLF5250-2</b>	5	6	250	3	10
	<b>DOLF5350-2</b>	5	6	350	3	10
	<b>DOPL4250-2</b>	4	8	250	7	10
	<b>DOLP4250-2</b>	4	8	250	8	10
	<b>DOLP4350-2</b>	4	8	350	8	10
	<b>DOLP5250-2</b>	5	6	250	8	10
	<b>DOLP5350-2</b>	5	6	350	8	10
	<b>PICK352L-2</b>	3.5	10	230	18	10
	<b>PICK352A-2</b>	4.0	8	300	8	10

## Vitruvian Handle



B 89ASH-LL

Luer Lock Hub

B 89TAH-60CC

60cc Toomey Syringe Hu



Luer Lock Aspiration Handle

Toomey Aspiration "Delrin" Handle

## Tip Styles



## Vitruvian Liposuction Cannula Popular Patterns

### Mercedes

B 8MER-215L	2 x 150mm	Vitruvian Handle
B 8MER-226L	2 x 260mm	Vitruvian Handle
B 8MER-215-SYR	2 x 150mm	Toomey Hub
B 8MER-326L	3 x 260mm	Vitruvian Handle
B 8MER-326-SYR	3 x 260mm	Toomey Hub
B 8MER-332L	3 x 320mm	Vitruvian Handle
B 8MER-415L	4 x 150mm	Vitruvian Handle
B 8MER-426-SYR	4 x 260mm	Toomey Hub
B 8MER-432L	4 x 320mm	Vitruvian Handle
B 8MER-526L	5 x 260mm	Vitruvian Handle

### Accelerator

B 8III-226L	2 x 260mm	Vitruvian Handle
B 8III-215-SYR	2 x 150mm	Toomey Hub
B 8III-315L	3 x 150mm	Vitruvian Handle
B 8III-326L	3 x 260mm	Vitruvian Handle
B 8III-332L	3 x 320mm	Vitruvian Handle
B 8III-426-SYR	4 x 260mm	Toomey Hub
B 8III-432L	4 x 320mm	Vitruvian Handle

### Fat Disruptor

B 8MBAT-315L	3 x 150mm	Vitruvian Handle
B 8MBAT-326-SYR	3 x 260mm	Toomey Hub

### Cobra

B 8COB-326L	3 x 260mm	Vitruvian Handle
B 8COB-426L	4 x 260mm	Vitruvian Handle

### Las Vegas

B 8LV-326-SYR	3 x 260mm	Toomey Hub
B 8LV-415L	4 x 150mm	Vitruvian Handle

### Spatula

B 8SPA-315L	3 x 150mm	Vitruvian Handle
B 8SPA-315-SYR	3 x 150mm	Toomey Hub
B 8SPA-415L	4 x 150mm	Vitruvian Handle
B 8SPA-415-SYR	4 x 150mm	Toomey Hub


### MicroAire Compatible Reusable Cannulae

Black & Black Surgical has developed a superior, reusable cannula, fabricated from the highest quality materials, using tight tolerances to minimise vibration.


Now available in a full line of sizes and patterns to meet all your liposuction needs.

Reverse switch for flow direction preference

- High quality materials and craftsmanship
- Made in the USA: U.S. Patent D950,045 / D950,046
- Cannula tubing connector designed for non-slip fit, for both aspiration and infiltration tubing
- Reusable / Autoclavable
- Validated cleaning and sterilization instructions



## Cannula Patterns



### Accelerator

B8B-III-326 PWR	3MM X 26CM STRAIGHT
B8B-III-426 PWR	4MM X 26CM STRAIGHT
B8B-III-432 PWR	4MM X 32CM STRAIGHT
B8B-III-436 PWR	4MM X 36CM STRAIGHT
B8B-III-526 PWR	5MM X 26CM STRAIGHT
B8B-III-532 PWR	5MM X 32CM STRAIGHT

### Mirrored Accelerator

B8B-MIII-326 PWR	3MM X 26CM STRAIGHT
B8B-MIII-332 PWR	3MM X 32CM STRAIGHT
B8B-MIII-426 PWR	4MM X 26CM STRAIGHT
B8B-MIII-432 PWR	4MM X 32CM STRAIGHT
B8B-MIII-526 PWR	5MM X 26CM STRAIGHT
B8B-MIII-532 PWR	5MM X 32CM STRAIGHT

### Basket

B8B-BAS-326 PWR	3MM X 26CM STRAIGHT
B8B-BAS-332 PWR	3MM X 32CM STRAIGHT
B8B-BAS-426 PWR	4MM X 26CM STRAIGHT
B8B-BAS-426 ANG PWR	4MM X 26CM ANGLED
B8B-BAS-432 PWR	4MM X 32CM STRAIGHT
B8B-BAS-432 ANG PWR	4MM X 32CM ANGLED
B8B-BAS-440 PWR	4MM X 40CM STRAIGHT
B8B-BAS-440 ANG PWR	4MM X 40CM ANGLED
B8B-BAS-526 PWR	5MM X 26CM STRAIGHT
B8B-BAS-526 ANG PWR	5MM X 26CM ANGLED
B8B-BAS-532 PWR	5MM X 32CM STRAIGHT
B8B-BAS-532 ANG PWR	5MM X 32CM ANGLED
B8B-BAS-536 PWR	5MM X 36CM STRAIGHT
B8B-BAS-536 ANG PWR	5MM X 36CM ANGLED

### Double Basket

B8B-DBAS-426 PWR	4MM X 26CM STRAIGHT
B8B-DBAS-432 PWR	4MM X 32CM STRAIGHT
B8B-DBAS-432 ANG PWR	4MM X 32CM ANGLED
B8B-DBAS-440 PWR	4MM X 40CM STRAIGHT
B8B-DBAS-440 ANG PWR	4MM X 40CM ANGLED
B8B-DBAS-526 PWR	5MM X 26CM STRAIGHT
B8B-DBAS-532 PWR	5MM X 32CM STRAIGHT

### Candy Cane

B8B-CC-326 PWR	3MM X 26CM STRAIGHT
B8B-CC-426 PWR	4MM X 26CM STRAIGHT
B8B-CC-432 PWR	4MM X 32CM STRAIGHT
B8B-CC-436 PWR	4MM X 36CM STRAIGHT
B8B-CC-526 PWR	5MM X 26CM STRAIGHT
B8B-CC-532 PWR	5MM X 32CM STRAIGHT

### Mercedes

B8B-III-326 PWR	3MM X 26CM STRAIGHT
B8B-III-426 PWR	4MM X 26CM STRAIGHT
B8B-III-432 PWR	4MM X 32CM STRAIGHT
B8B-III-436 PWR	4MM X 36CM STRAIGHT
B8B-III-526 PWR	5MM X 26CM STRAIGHT
B8B-III-532 PWR	5MM X 32CM STRAIGHT

### Double Mercedes

B8B-DMER-326 PWR	3MM X 26CM STRAIGHT
B8B-DMER-426 PWR	4MM X 26CM STRAIGHT
B8B-DMER-432 PWR	4MM X 32CM STRAIGHT
B8B-DMER-436 PWR	4MM X 36CM STRAIGHT
B8B-DMER-526 PWR	5MM X 26CM STRAIGHT
B8B-DMER-532 PWR	5MM X 32CM STRAIGHT
B8B-DMER-536 PWR	5MM X 36CM STRAIGHT

### Helix 3-Port

B8B-HLX-326 PWR	3MM X 26CM STRAIGHT
B8B-HLX-332 PWR	3MM X 32CM STRAIGHT
B8B-HLX-426 PWR	4MM X 26CM STRAIGHT
B8B-HLX-432 PWR	4MM X 32CM STRAIGHT
B8B-HLX-436 PWR	4MM X 36CM STRAIGHT
B8B-HLX-526 PWR	5MM X 26CM STRAIGHT
B8B-HLX-532 PWR	5MM X 32CM STRAIGHT
B8B-HLX-536 PWR	5MM X 36CM STRAIGHT

### Large Volume Harvester


B8B-LVH-326 PWR	3MM X 26CM STRAIGHT
B8B-LVH-332 PWR	3MM X 32CM STRAIGHT
B8B-LVH-426 PWR	4MM X 26CM STRAIGHT
B8B-LVH-432 PWR	4MM X 32CM STRAIGHT
B8B-LVH-440 PWR	4MM X 40CM STRAIGHT
B8B-LVH-526 PWR	5MM X 26CM STRAIGHT
B8B-LVH-532 PWR	5MM X 32CM STRAIGHT

### Spatula

B8B-SPA-315 PWR	3MM X 15CM STRAIGHT
B8B-SPA-415 PWR	4MM X 15CM STRAIGHT


## Single Use Fat Harvesting Cannula

Cannulas specifically designed to perform an accurate fat SIZING during harvesting thanks to our patented diamond-shaped ports with 360° cutting edges.



- MEGAFAT
- Small volume
- LL syringe
- Suction 360°

- MACROFAT
- Med/Large volume
- LL syringe / Pump
- Suction precision






















Metallic Luer/Lock

A range of cannula specifically designed to perform accurate fat SIZING during harvesting due to the diamond shaped cutting ports with 360° cutting edges, when selecting the appropriate negative pressure to harvest.

	REF	Ø mm	G	↔ mm	Holes	Box
<b>HARVESTING</b>						
	<b>PICK3518-2</b>	3.5	10	180	5	10
	<b>PICK352E-2</b>	2.0	14	150	8	10
	<b>PICK352L-2</b>	3.5	10	230	18	10
Luer lock with handle	<b>PICK352A-2</b>	4.0	8	300	8	10
	<b>CP3015</b>	3	11	150	2	10
	<b>CP3023</b>	3	8	230	2	10
	<b>CP161053HS</b>	1.65	16	105	3	10
	<b>CP141253HS</b>	2.1	14	125	3	10
	<b>CP30153S</b>	3	11	150	3	10
	<b>CP31503M</b>	3	11	150	3	10
	<b>CP32303M</b>	3	11	230	3	10
	<b>CP42303M</b>	4	8	230	3	10
	<b>CP32303T</b>	3	11	230	3	10
	<b>CP42303T</b>	4	8	230	3	10
	<b>CP30154H</b>	3	11	230	4	10
	<b>CP141256HMP</b>	2.1	14	125	6	10
	<b>CP141506HS</b>	2.1	14	150	6	10
	<b>CP141508H</b>	2.1	14	150	8	10
	<b>CP113102 (Rough)</b>	3	11	300	12	10
	<b>CP301512R (Rough)</b>	3	11	150	18	10
	<b>CP301524H</b>	3	11	150	24	10
	<b>CP302324H</b>	3	11	230	24	10
	<b>CP122007H</b>	2.8	12	200	7	10
	<b>CP123007H</b>	2.8	12	300	7	10
	<b>DOTO 0155-2 HANDLE</b>					10

## Single Use Fat Injection Needles

	REF	Ø mm	G	↔ mm	Holes	Box
<b>FAT INJECTION</b>						
	<b>FIL1270S-2</b>	3.5	10	250	1	10
	<b>FILU2515-2</b>	2.5	12	200	1	10
	<b>FIL2590S-2</b>	2.7	12	200	1	10
	<b>F9027150-2</b>	2.7	12	150	1	10
	<b>FIL25C5S-2</b>	2.0	14	150	1	10
	<b>FIL1670S-2</b>	2.0	14	150	1	10
	<b>F9020100-2</b>	2.0	14	100	1	10
	<b>F9015100-2</b>	1.5	17	100	1	10
	<b>F9015070-2</b>	1.5	17	70	1	10
	<b>F9012100-2</b>	1.2	18	100	1	10
	<b>F9012070-2</b>	1.2	18	70	1	10
	<b>CF151045</b>	1.65	16	100	1	10
	<b>CF201545</b>	2.10	14	150	1	10
	<b>CF202345</b>	2.10	14	230	1	10
	<b>CF10790</b>	1.05	19	70	1	10
	<b>CF15790</b>	1.65	16	70	1	10
	<b>CF16890</b>	1.65	16	80	1	10
	<b>CF151090</b>	1.65	16	100	1	10
	<b>CF201090</b>	2.10	14	100	1	10
	<b>CF201590</b>	2.10	14	150	1	10
	<b>CF202290</b>	2.10	14	220	1	10
	<b>CF301590</b>	3.00	11	150	1	10
	<b>CP1690V</b>	1.65	16	90	1	10
	<b>CP1490V</b>	2.10	14	90	1	10
 Curved	<b>CP142090CSP</b>	2	14	200	1	10

	REF	Ø mm	G	↔ mm	Holes	Box
<b>FAT INJECTION</b>						
	<b>MI18110KR</b>	1.2	18	110	1	20
	<b>MI1890KR</b>	1.2	18	90	1	20
	<b>MI1870KR</b>	1.2	18	70	1	20
	<b>F9010100</b>	1	19	100	1	20
	<b>F9010070</b>	1	19	70	1	20
	<b>MI2170KR</b>	0.81	21	70	1	20
	<b>MI2140KR</b>	0.81	21	40	1	20
	<b>F9008060</b>	0.80	21	60	1	20
	<b>MI2290KR</b>	0.70	22	90	1	20
	<b>MI2270KR</b>	0.70	22	70	1	20
	<b>MI2250KR</b>	0.70	22	50	1	20
	<b>MI2370KR</b>	0.63	23	70	1	20
	<b>MI2350KR</b>	0.63	23	50	1	20
	<b>MI2450KR</b>	0.55	24	50	1	20
	<b>MI2438KR</b>	0.55	24	38	1	20
	<b>MI2550KR</b>	0.50	25	50	1	20
	<b>MI2650KR</b>	0.45	26	50	1	20
	<b>MI2638KR</b>	0.45	26	38	1	20
	<b>MI2625KR</b>	0.45	26	25	1	20
	<b>MI2740KR</b>	0.40	27	40	1	20






**Centrifuge for Fat Grafting**

CES-10 variable speed, time adjustable medical centrifuge for use with fat harvesting preparation and Platelet Rich Plasma. Can be used in the doctors office or in the operating theatre. The new CES-10 comes with a locking lid and 8 metal syringe holders / tubes - that are designed to hold 8 x 10ml BD luer lock syringes - with plunger REMOVED. The 8 metal syringe tubes/holders are autoclavable.

Mains powered and fitted with noise and vibration reducing pads.

## Reusable Fat Injection

	REF	∅ mm	↔ mm	Shape
<b>FAT INJECTION</b>				
 <p>Style I</p>	<b>B8FIC-I-1415-LL</b>	14	150	Straight
	<b>B8FIC-I-1420-LL</b>	14	200	Straight
	<b>B8FIC-I-1425-LL</b>	14	250	Straight
	<b>B8FIC-I-167-LL</b>	16	70	Straight
	<b>B8FIC-I-169-LL</b>	16	90	Straight
	<b>B8FIC-I-187-LL</b>	18	70	Straight
	<b>B8FIC-I-189-LL</b>	18	90	Straight
	<b>B8FIC-I-195-LL</b>	19	50	Straight
	<b>B8FIC-I-CC1415-LL</b>	14	150	Concave
	<b>B8FIC-I-CC1610-LL</b>	16	100	Concave
	<b>B8FIC-I-CC169-LL</b>	16	90	Concave
	<b>B8FIC-I-CV1415-LL</b>	14	150	Convex
	<b>B8FIC-I-CV1610-LL</b>	16	100	Convex
	<b>B8FIC-I-CV169-LL</b>	16	90	Convex
 <p>Style II</p>	<b>B8FIC-II-167-LL</b>	16	70	Straight
	<b>B8FIC-II-169-LL</b>	16	90	Straight
	<b>B8FIC-II-187-LL</b>	18	70	Straight
	<b>B8FIC-II-189-LL</b>	18	90	Straight
	<b>B8FIC-II-CC1610-LL</b>	16	100	Concave
	<b>B8FIC-II-CC169-LL</b>	16	90	Concave
	<b>B8FIC-II-CV1610-LL</b>	16	100	Convex
	<b>B8FIC-II-CV169-LL</b>	16	90	Convex
 <p>'V' Style</p>	<b>B8FIC-V-1215-LL</b>	12	150	Straight
	<b>B8FIC-V-1225-LL</b>	12	250	Straight
	<b>B8FIC-V-1415-LL</b>	14	150	Straight
	<b>B8FIC-V-1425-LL</b>	14	250	Straight

## Small Volume Lipofilling


st'rim®

MICROFAT grafting kit - REF PICS2130-2



SIZING



360° SUCTION



MEDIUM FLOW



	TIP	A	B1	B2	B3	C1	C2	REF	Box
St'rim	∅ ↔mm	2.0 mm 130	19G 50	21G 60	21G 40	14G 40	21G 40	PICS2130-2	5

## Lipocube




**Nanocube** "All in one" Fat Grafting - 4 ports for small volume fat grafting procedures;

Sterile, single-use device for re-sizing harvested fat particles into mini, micro and nanofat grafts for specific indications;

Supplied sterile, pack size each (1) REF L-NANO

### INDICATIONS




Milifat (port 1-2) 2.4mm<

Deep fat compartments, Fat-Lifting, Bone onlay graft



Microfat (port 3) 1.2mm<

Superficial fat compartments




Nanofat (port 4) 600 microns <

Subdermal, Intradermal, Masotherapy, Micro-needling, topical

st'rim®

**MACROFAT** grafting kit - **REF PICS3523-2**



	TIP	A	B1	B2	B3	C1	C2	REF	Box
Macro St'rim	∅ ←mm	3.5 mm 220	12G 150	14G 150	17G 100	17G 70		PICS3523-2	5

**DOTO 0155-2 HANDLE**

To convert Strim Macro kit harvesting cannula from luer lock to liposuction handle - Box of 10



# NLF - Fat Collection System

Sterile, single use kit for fat collection during liposuction with double filter to allow swift fluid removal. There is an easy access non-return, luer lock valve for collection of high quality fat cells.

**REF NLF2SYSTEM, NLF2STAND**

- Completely closed and sterile circuit
- Biocompatible
- Single-use
- Compatible with the manual liposuction cannulae and power assisted devices
- Easy to use
- In compliance with the European standards - CE marking / Medical device class IIa



## Accessories



### Toomey Syringe Luer Lock Adaptor - Reuseable

Convert the Toomey syringe with this reuseable, autoclavable adaptor to connect Luer Lock Fat Harvesting Cannula.

Order Ref; B89TFL



### Toomey Syringe

Sterile 50/60ml Toomey style syringe to fit with Toomey style cannula. Available in a box of 5 syringes. Not CE Marked.

Order Ref; 89-SYR-60TT Monoject syringe



### Syringe Snap Lock

Locks to help maintain syringe pressure, these “snap” locks connect to the back finger rest of the syringe case. Lock is maintained by twisting the plunger of the syringe to hold in place. Reuseable and autoclavable, these are available for 10ml and 60ml syringes and supplied individually.

Order Ref; B89TL-10CC 10ml Syringe Lock

Order Ref; B89TL-60CC 60ml Syringe Lock



### Re-useable Syringe Snap

Order Ref; B89SL-60CC 60 ml

Order Ref; B89SL-10CC 10 ml



### Syringe Racks Reuseable

Reuseable and autoclavable steel syringe racks for holding syringes of harvested fat prior to injection. Available in the following sizes:

Order Ref;	Rack	Syringe Size Holes
POR6383	3 x 6	Syringes 1 - 10ml
POR36113	3 x 6	Syringes 1 - 50ml



### Luer Lock Syringe Cap

Fits all standard luer lock syringes.

Order Ref; MC1711A (box of 100)



### Luer Lock Transfer Adaptor

Male to male luer lock syringe connector (2.4mm) to transfer fat from one syringe to another.

Order Ref; LLF24-100 (box of 100)

Order Ref; B89LL3 (Reusable Stainless Steel)




## Female Garments




### 9" Abdominal Binder

Using a unique cotton fabric, 9» (23cm) panel structure, with a full length Velcro fastening, the binder can be fitted to all sizes and the amount of compression adjusted easily and on a daily base, that will allow you to remain in comfort whilst you assist and recover from surgery. - [Order Ref; 355-9](#)




### Full body Girdle - Zippered

Stage 1 Zippered Garments allow greater ease when getting the garment on / off particularly important in the early stages of recovery after surgery. The Stage 1 garments are designed to allow for greater expansion to accommodate post-surgical swelling with appropriate compression to match. - [Order Ref; 855](#)



### Abdominal Girdle - Zippered

The Female Abdominal Girdle Compression Cosmetic Surgery Garment (with Zipper) - Stage One features zippers for easy application, adjustable shoulder straps, targeted compression for the abdominal area, reinforced flat seams, gripper lace to prevent roll up, and a reinforced open crotch. - [Order Ref; 850](#)




### Above-knee Girdle - Zippered

The Female Above Knee Girdle Compression Garment (with Zipper) - Stage One Garment improves blood circulation, minimises swelling after the procedure, flushes the body out of potentially harmful fluids, accelerates the healing process, and allows the patient to return to daily routines sooner. - [Order Ref; 753](#)



### Abdominal Brief - Zippered

The Female Abdominal Brief Cosmetic Surgery Compression Garment (with Zipper)- Stage One Garment improves blood circulation, minimizes swelling after the procedure, flushes the body out of potentially harmful fluids, accelerates the healing process, and allows the patient to return to daily routines sooner. - [Order Ref; 75](#)



### Compression Arm Sleeve With Adjustable Cotton Knit Bra

Made with 3/4 inch length compression arm sleeves, built-in seamless cotton knit bra, convenient front and back hook and eye closure, and no wiring and no seams at the surgical areas for maximum comfort and to avoid irritation, the Plastic Surgery Compression Arm Sleeves w/Adjustable Cotton Knit Bra - Moderate Support is ideal for a patient having multiple procedures. - [Order Ref; B943](#)

### Male Garments



#### 12" Triple Panel Binder

This 12" triple panel binder is perfect following abdominal and flank procedures. It is made from strong, durable stretch fabric which is lightweight, cool and comfortable. - [Order Ref; 340-12](#)



#### Zippered Abdominal / Chest Garment

The Zippered Abdominal/Chest Garment is ideal for chest, abdominal & flank procedures. The breathable powernet fabric retains compressive memory is washable and fast drying and is lightweight, cool and comfortable. - [Order Ref; V840](#)



#### Zippered Above-knee Abdominal Garment

This garment is ideal following abdominal, flank, hip & thigh procedures. It has been carefully designed to incorporate side zip with cotton guard tape & high waist design with 3X wide waistband for ultimate comfort post operation. - [Order Ref; 743](#)



#### Adjustable Compression Vest

The Male Cosmetic Surgery Compression Vest provides compression in the abdomen, chest, and flank areas and is appropriate for liposuction, gynecomastia, and pectoral implants. This garment improves blood circulation, minimizes swelling after the procedure, flushes the body out of potentially harmful fluids, accelerates the healing process, and allows the patient to return to daily routines sooner. - [Order Ref; 640](#)

### Accessories



#### Breast Icing Cryotherapy Wrap

The Breast Icing wrap is a cryotherapy kit with two gel cold packs for post surgery recovery. The packs are frozen for two hours, worn for 20 minutes, and can be refrozen for repeated use. [Ref: EST60-SEIN](#)



#### Epifoam Compression Pads

Epifoam offers patients an important evolution in liposuction recovery. Providing a superior platform over which to apply post-operative compression garments. - [Order Ref; EPF30](#)

## Our Garments Range

### Female Garments

#### Abdominal & Mid Body

- 9" Abdominal Binder - [Order Ref; 355-9](#)
- 6" Single Panel Binder - [Order Ref; 350-6](#)
- 9" Double Panel Binder - [Order Ref; 350-9](#)
- Torso/Brief Zippered - Stage 1 - [Order Ref; 85](#)
- Torso/Brief Non-Zip - Stage 2 - [Order Ref; 165](#)
- High back Torso/Brief With Bra Zippered - Stage 1 - [Order Ref; B85](#)
- High back Torso/Brief With Bra Non-Zip - Stage 2 - [Order Ref; B165](#)

#### Lower Body

- Labia Support Thong - [Order Ref; 1052](#)
- Labia Support Brief - [Order Ref; 1054](#)
- Below-knee Girdle Zippered - Stage 1 - [Order Ref; 755](#)
- Below-knee Girdle Non-Zip - Stage 2 - [Order Ref; 1455](#)
- Above-knee Girdle Zippered - Stage 1 - [Order Ref; 753](#)
- Above-knee Girdle Non-Zip - Stage 2 - [Order Ref; 1453](#)
- Tram Flap Girdle Zippered - Stage 1 - [Order Ref; 750](#)
- Tram Flap Girdle Non-Zip - Stage 2 - [Order Ref; 1450](#)

### Female Garments

#### Abdominal

- 9" Abdominal Binder - [Order Ref; 345-9](#)
- 12" Triple Panel Binder - [Order Ref; 340-12](#)
- 9" Double Panel Binder - [Order Ref; 340-9](#)

#### Upper Body

- Dorsocervical Male Vest - [Order Ref; DC-64](#)
- Zippered Compression Vest With Arms - Stage 1 - [Order Ref; 642](#)
- Non-zip Compression Tank Top - Stage 2 - [Order Ref; 1241](#)
- Zippered Compression Tank Top - Stage 1 - [Order Ref; 641](#)
- Compression Vest - Stage 2 - [Order Ref; 1240](#)
- Adjustable Compression Vest - Stage 1 - [Order Ref; 640](#)

#### Upper Body

- Dorsocervical Female Garment - [Order Ref; DC-65](#)
- Long Arm Sleeve With Lace - [Order Ref; L933](#)
- Compression Arm Sleeve With Adjustable Bodice - [Order Ref; 943](#)
- Compression Arm Sleeve With Adjustable Bra - [Order Ref; B943](#)
- Long Arm Sleeve Bolero - [Order Ref; BV933](#)
- Short Arm Sleeve Bolero - [Order Ref; BV930](#)

#### Upper Body




- Full Body Girdle Zippered - Stage 1 - [Order Ref; 855](#)
- Full Body Girdle Non-Zip - Stage 2 - [Order Ref; 1655](#)
- High-Back Body Girdle Zippered - Stage 1 - [Order Ref; 853-H](#)
- High-Back Body Girdle Non-Zip - Stage 2 - [Order Ref; 1653-H](#)
- Below-knee Body Girdle Zippered - Stage 1 - [Order Ref; BA855](#)
- Below-knee Body Girdle Non-Zip - Stage 2 - [Order Ref; B1655](#)
- Above-knee Body Girdle Zippered - Stage 1 - [Order Ref; B853](#)
- Above-knee Body Girdle Non-Zip - Stage 2 - [Order Ref; B1653](#)

#### Lower Body

- Scrotum Support Brief - [Order Ref; 1045](#)
- Zippered Above-knee Abdominal Garment - Stage 1 - [Order Ref; 753](#)
- Zippered Gluteus Garment - Stage 1 - [Order Ref; G740](#)

#### Lower Body

- Non-zip Abdominal Garment - Stage 2 - [Order Ref; 1640](#)
- Zippered Abdominal Garment - Stage 1 - [Order Ref; 840](#)
- Non-zip Abdominal/Chest Garment - Stage 2 - [Order Ref; V1640](#)
- Zippered Abdominal/Chest Garment - Stage 1 - [Order Ref; V840](#)
- Zippered Bodysuit - Stage 1 - [Order Ref; BS843](#)





For post surgical support garments, visit  
[www.esshop.co.uk](http://www.esshop.co.uk)



**NLF2 kit**



**Lipocube Nanocube**

**Eurosurgical Limited**  
Merrow Business Centre, Guildford,  
Surrey GU4 7WA, United Kingdom

+44 (0) 1483 456 007  
[sales@eurosurgical.co.uk](mailto:sales@eurosurgical.co.uk)  
[www.eurosurgical.co.uk](http://www.eurosurgical.co.uk)

V003 01 26